

NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device.


The SSCP is not intended to replace the Instructions for Use as the main document to ensure the safe use of the device, nor is it intended to provide diagnostic or therapeutic suggestions to intended users or patients.

The following information is intended for users/healthcare professionals.

1. Device identification and general information	
Device trade name(s)	<u>NuMED Stent Placement Family</u> BIB® (Balloon in Balloon) Stent Placement Catheter
Model Number	<u>NuMED Stent Placement Family – Model 1500</u> BIB – Model 420.1
Manufacturer's name and address	NuMED, Inc. 2880 Main Street Hopkinton, NY 12965 USA
Manufacturer's single registration number (SRN)	US-MF-000010948
Basic UDI-DI	08877141500SV
Medical device nomenclature description / text	EMDN – C010402020102 – Cardiocirculatory System Devices, Stent Positioning Vascular Balloon Dilatation Catheters
Class of device	III
Year when first certificate (CE) was issued	2003
Authorised Representative (AR)	EVOMED, S.L.U. Ctra. Torrejón-Ajalvir Km. 5,2 28864 Ajalvir (Madrid), Spain
AR SRN	ES-AR-000047116
Notified Body	SGS Belgium NV
Notified Body ID number	1639

2. Intended use of the device	
Indications for use	Indicated for stent placement in vessels over 8mm in diameter.
Contraindications and/or limitations	There are no contraindications to the use of the BIB Stent Placement Catheters. Limitations – only for use in vessels over 8mm in diameter.

3. Device description	
Description of the device	The BIB Stent Placement Catheter is triaxial in construction with two lumens being used to inflate the balloon while one lumen is being used for tracking over a guidewire. The inner balloon is 1/2 of the outer balloon diameter and 1 cm shorter. The purpose of the double balloon catheter is to apply an incremental inflation for the purpose of dilating a stent. The inner balloon provides initial expansion of the stent and also acts as a tool to hold the stent on the catheter prior to the outer balloon being inflated.



NuMED


Summary of Safety and Clinical Performance

SSCP – Stent Placement

	<p>The outer balloon is then inflated, providing the remainder of the expansion. There are radiopaque platinum marker bands under the balloon shoulders, to aid during placement.</p> <p>The balloons are designed to inflate to the diameter and length listed on the label at a specific pressure. Thus, it is recommended that the device be used in conjunction with a mechanism to monitor pressure, an inflation device with pressure gauge. This device is also designed to be used with an appropriately sized introducer and guidewire, and stent.</p> <p>Both the inner and outer balloon size is $\pm 10\%$ at the Rated Burst Pressure (RBP) and the RBP is not to be exceeded.</p> <p>The catheters are supplied sterile, by ethylene oxide gas, and are intended for single use only. The catheters are invasive and intended for transient use (continuous use of <60 minutes) on patients.</p>
Reference to previous generation(s) or variants	The BIB also comes in variants with a stent pre-mounted on it: Mounted CP Stent and Covered Mounted CP Stent. These other variants are covered by NuMED's CoA Stent and RVOT/CoA SSCPs.
Accessories which are intended to be used in combination with the device	Mandrel – for use when mounting the stent on the catheter.
Description of any other devices and products which are intended to be used in combination with the device	This device is designed to be used with a stent, guidewire, introducer sheath, inflation device with pressure gauge, and balloon inflation media.

4. Risks and Warning

Residual risks and undesirable effects	<p>Side-effects reported in the literature are stent flaring and stent migration.</p> <p>All risks identified in the clinical literature as well as the risks detected from the Post Market Surveillance or from clinical data generated and held by the Sponsor, have been considered by the risk management process.</p> <p>All significant risks were considered, mitigated as far as possible (AFAP), and are acceptable in regard to the clinical benefit of the device.</p> <p>Identified clinical residual risks/undesirable side-effects for the Stent Placement Catheters are:</p> <p>POTENTIAL COMPLICATIONS/ADVERSE EFFECTS</p> <p>NOTE: Circumferential tear of the delivery balloon catheter prior to complete expansion of the stent may cause the balloon to become tethered to the stent, requiring surgical removal. In case of rupture of an adequately sized balloon after stent expansion, it can be withdrawn and a new balloon catheter exchanged over a guidewire to complete expansion of the stent.</p> <p>Cardiac catheterization carries certain risks. Potential complications & adverse effects associated with device use and indication include:</p> <ul style="list-style-type: none"> • Stent misplacement • Stent migration • Minor hematoma • Intraluminal thrombosis • Pseudoaneurysm • AV fistula formation • Bleeding • Sepsis/infection • Distal thromboemboli • Death
--	--




NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

	<ul style="list-style-type: none"> • Vessel rupture • Cerebrovascular incident • Hematoma requiring repair • Femoral Artery Injury
Warning and Precautions	<p>The following Warnings and Precautions have been identified and are called out in the Instructions for Use:</p> <p>WARNINGS</p> <ul style="list-style-type: none"> • Do not exceed the RBP. An inflation device with pressure gauge is recommended to monitor pressure. Pressure in excess of the RBP can cause balloon rupture and potential inability to withdraw the catheter through the introducer sheath. • Confirm that the distal end of the introducing sheath is at least 2.5cm back from the most proximal image markers before inflating the outer balloon. Failure to do so may stretch the outer tubing and severely hinder balloon deflation. • Use two appropriate size inflation devices with pressure gauges for inflation. • Do not advance the guidewire, balloon dilatation catheter, or any other component if resistance is met, without first determining the cause and taking remedial action. • When the stent is crimped onto a balloon delivery catheter, the maximum balloon inflation pressure must not exceed the recommended inflation pressure specified in the manufacturer's instructions. • This catheter is not recommended for pressure measurement or fluid injection. • Do not remove the guidewire from the catheter at any time during the procedure except when the procedure has been completed. • This device is intended for single use only. Do not resterilize and/or reuse it, as this can potentially result in compromised device performance and increased risk of cross-contamination. <p>PRECAUTIONS</p> <ul style="list-style-type: none"> • The BIB Stent Placement balloon catheter was tested with the NuMED Cheatham Platinum (CP) Stent (bare & covered). • Dilatation procedures should be conducted under fluoroscopic guidance with appropriate x-ray equipment. • Stents are delicate devices. Care should be exercised while handling to help prevent the possibility of breakage. • Careful attention must be paid to the maintenance of tight catheter connections by aspiration before proceeding to avoid air introduction into the system. • The inflation diameter of the balloon used during stent delivery should approximate the diameter of the obstructive vessel and the intended implant site. • Under no circumstances should any portion of the catheter system be advanced against resistance. The cause of the resistance should be identified with fluoroscopy and action taken to remedy the problem. • If resistance is felt upon removal, then the balloon, guidewire and the sheath should be removed together as a unit, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the balloon catheter and sheath as a unit and withdrawing both together, using a gentle twisting motion combined with traction. • Before removing the catheter from the sheath, it is very important that the balloon is completely deflated. • Proper functioning of the catheter depends on its integrity. Care should be used when handling the catheter. Damage may result from kinking, stretching, or forceful wiping of the catheter.
Other relevant aspects of safety, including a summary of any field safety corrective actions (FSCA including FSN) if applicable	<p>There have not been any Field Safety Corrective Actions or Field Safety Notices for the BIB Stent Placement Catheter.</p>



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

5. Summary of clinical evaluation and post-market clinical follow-up (PMCF)


Summary of clinical data related to equivalent device:

An equivalent device was not used for the clinical evaluation.

Summary of clinical data from conducted investigations of the device:

The COAST study was specifically initiated to investigate clinical efficacy and safety of the CP Stent, with the BIB Stent Placement Catheter. The result of this study demonstrates that a double balloon significantly reduces known risks associated with single balloon stent placement procedures. Devices used in the studies were provided both in unmounted configurations.

1. Study name: COAST										
Appraisal										
Level of Evidence	Study Method/Design		Question Applied			Oxford LOE 2011				
	Prospective, multicenter, single arm interventional, open label study.		To evaluate the intermediate results of CP Stent (Bare and Covered) to treat native and recurrent CoA in selected children, adolescents and adult.			1	2	3	4	5
Suitability	Relevant Data					Grading				
Device	- CP Stents (Bare and Covered) and BIB - BIB was manually crimped on the Bare Stents - Covered Stents were provided pre-mounted on the BIB					D1	D2		D3	
Application	- CoA (native and recurrent)					A1	A2		A3	
Patient	- Patients with native or recurrent CoA - Sampling: n=105 - Mean age: 16 (range: 8 to 52) years old - Sex: 73M; 32F					P1	P2		P3	
Report	- High quality					R1	R2		R3	
Suitability Grade (Range 4-12)						4				
Data Contribution	Relevant Data					Grading				
Outcomes/Endpoints	- Blood pressure gradient - Coarctation minimum diameter: cardiac catheterization before and after CP Stent placement - Safety					Yes 1		No 2		
Follow-up	- 12 - 24 months with certain assessments extending to 5 years.					Yes 1		No 2		
Statistical analysis	- Descriptive statistics are presented as mean ± SD or median (minimum–maximum). Bivariate comparisons of preimplantation and post implantation catheterization data and subsequent blood pressures were performed with the paired t test. Comparison of means or proportions between populations were performed by unpaired t test or Wilcoxon rank-sum test based on distribution and the Fisher exact test, respectively. Multivariable analysis of dichotomous outcome variables was performed with logistic regression. Analysis of time-dependent occurrences was presented graphically with Kaplan-Meier plots and analyzed statistically by the log-rank test. Predictors of time-dependent outcomes such as reintervention were obtained from Cox proportional hazards modeling.					Yes 1		No 2		
Clinical significance	- The CP stent is safe and associated with persistent relief of aortic obstruction. Stent fracture and progression of fracture occur but have not resulted in clinically important sequelae. Reintervention is common and related to early and late aortic wall injury and need for re-expansion of small-diameter stents.					Yes 1		No 2		
Data Contribution Grade (Range 4-8)						4				
Overall S&P Appraisal, Disposition and Weighting										
S&P Grade (Range 9-25)	LOE (2) + Suitability (4) + Data Contribution (4) = 10		Disposition and Weighting (select)		Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21 Excluded, 22-25					
Purpose: to provide information that will support labeling of both the CP Stent (Bare and Covered) to treat native and recurrent CoA in selected children, adolescents and adult										
Clinical Study Methodology: Single arm interventional study (open label). The COAST is a prospective, multicenter, single-arm clinical study										



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

involving 19 pediatric cardiology centers in the United States. The study includes patients with native or recurrent CoA treated by physicians at the participating institutions. A total of 105 patients underwent attempted implantation, with 104 successes.

Reference to the clinical study plan (and amendment) n°: NCT00552812

Investigation Site: 19 pediatric cardiology centers in United States Reference to Approved Consent Forms: N/A	Ethics Committee Approvals: Institutional Review Board approvals from all participating institutions Reference to Document n°: N/A	Regulatory Authority Approvals: Investigational Device Exemption from US FDA (August 3, 2007) Reference to Documents n°: G060057
--	--	--

Patient Population: Patients with native or recurrent CoA. A total of 105 patients underwent attempted implantation, median age 16 years (range from 8 to 52 years) and with 69.5% male.

Clinical Study Results: Results held by Sponsor.


Purpose	Criteria	Results
Performance	Blood pressure gradient and coarctation minimum diameter: cardiac catheterization before and after CP Stent placement	Average systolic blood pressure difference (mmHg) changed from 29±14 mmHg at baseline to -3±15 mmHg at 24 months follow-up. The Coarctation minimum diameter reported at 7.9 ± 2.7mm at baseline to 14±3 mm after implantation.
Safety	Adverse events	No serious adverse events reported, 7% of the patients experienced somewhat serious events. Aortic aneurysms (n=6): 5 were successfully treated with covered stent placement, and 1 resolved without intervention. Stent fractures were seen in 2 patients after one year, 11 patients at two years and 12 additional fractures above 2 years.

Reference to the Clinical Study Report n°: NCT00552812

Device Used: Bare CP Stent and BIB catheter; covered stents were available in case of aortic wall injury

Conclusion: The CP stent is safe and associated with persistent relief of aortic obstruction. Stent fracture and progression of fracture occur but have not resulted in clinically important sequelae. Reintervention is common and related to early and late aortic wall injury and need for re-expansion of small-diameter stents.

Clinical Publication: Meadows J, Minahan M, McElhinney DB, McEnaney K, Ringel R. Intermediate Outcomes in the Prospective, Multicenter Coarctation of the Aorta Stent Trial (COAST). Circulation 131, 1656-1664 (2015) (17)



NuMED


Summary of Safety and Clinical Performance

SSCP – Stent Placement

Summary of clinical data from other sources:

First Author (Year)	Appraisal/Results										
1. Sohrabi et al. (2014)	Safety & Performance Appraisal										
	Level of Evidence		Study Method/Design		Question Applied			Oxford LOE 2011			
			Prospective randomized controlled trial.		To evaluate outcomes of treatment with Covered versus Bare NuMED CP Stents.			1	2	3	4
	Suitability		Relevant Data						Grading		
	Device		- NuMED CP Stent (Bare and Covered) - BIB						D1	D2	D3
	Application		- Severe native CoA						A1	A2	A3
	Patient		- Patients with severe native CoA - Sampling: n=120 (60 CP Stents versus 60 CP Stents, Covered) - Mean age: 23.6±10.99 (range 12 to 58) years - Sex: 79 M; 41 F						P1	P2	P3
	Report		- High quality.						R1	R2	R3
	Suitability Grade (Range 4-12)							4			
	Data Contribution		Relevant Data						Grading		
	Outcomes/Endpoints		- Procedural success - Reduction in systolic blood pressure gradient - Reduction in mean diameter of coarctation segment - Adverse effects						Yes 1		No 2
Follow-up		- 31.1 ± 19.2 months						Yes 1		No 2	
Statistical analysis		- A p-value <0.05 was considered significant.						Yes 1		No 2	
Clinical significance		- Implanting CP Stent (Bare) and CP Stent (Covered) have very high success rates with remarkable hemodynamic effects in severe native CoA patients, with no significant complication during the procedure and hospitalization. - Patients undergoing CP Stent (Covered) implantation experienced a non-significantly lower re-coarctation rate and a higher occurrence of pseudoaneurysm formation with respect to CP Stent (Bare) stenting during follow-up. - In both groups, blood pressure was significantly reduced after intervention. - These findings indicate that CoA stenting is a safe procedure.						Yes 1		No 2	
Data Contribution Grade (Range 4-8)							4				
Overall S&P Appraisal, Disposition and Weighting											
S&P Grade (Range 9-25)		LOE (2) + Suitability (4) + Data Contribution (4) = 10			Disposition and Weighting (select)		Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21 Excluded, 22-25				

Contribution	
S&P	x
SOA	




NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

	Relevant S&P Results									
	Safety data	<ul style="list-style-type: none">- Pseudoaneurysms: 0 (CP Stent, Bare) versus 2 (CP Stent, Covered)- Mortality: 1 (CP Stent, Bare) versus 0 (CP Stent, Covered)								
	Performance data	<ul style="list-style-type: none">- Successful placement: successful in all patients- Mean systolic blood pressure gradient reduction: from 54.61 (CP Stent, Bare) and 54.42 (CP Stent, Covered) to 3.47 and 3.36 mmHg respectively; no significant difference between the two types of stents, P<0.001- Mean diameter of coarctation segment reduction: From 3.34 (CP Stent, Bare) and 3.30 (CP Stent, Covered) to 16.07 and 15.82 mm respectively; no significant difference between the two types of stents, P<0.001- Recurring coarctation: 4 (CP Stent, Bare) versus 0 (CP Stent, Covered), non-significant								
	Benefits/claims data	<ul style="list-style-type: none">- Reduction in mean systolic blood pressure gradient- Reduction in diameter of coarctation segment								
	Strengths	<ul style="list-style-type: none">- The CP Stent was hand-crimped down onto a balloon-in-balloon catheter (NuMED), which allows a precise and safe stent delivery								
	Weaknesses/ Potential bias	<ul style="list-style-type: none">- Although the first randomized clinical trial in this respect, study was limited in some aspects. First, during follow-up, patients did not undergo 24-hour ambulatory blood pressure monitoring, which could have diagnosed the normotensive state more accurately. Second, evaluation of the blood pressure response during exercise testing could have been more valuable in defining the procedure outcome.								
2. Vanagt et al. (2014)	Safety & Performance Appraisal									
	Level of Evidence	Study Method/Design	Question Applied			Oxford LOE 2011				
		Single-center retrospective study (CHD database of all CP Stent, Covered, during 2003-2012)	To evaluate possibilities and safety of CP Stent (Covered) in CHD.			1	2	3	4	5
	Suitability	Relevant Data				Grading				
	Device	<ul style="list-style-type: none">- CP Stent (Covered)- BIB				D1	D2	D3		
	Application	CoA and RVOT pre-stenting for percutaneous revaluation				A1	A2	A3		
	Patient	<ul style="list-style-type: none">- Patients with CoA and RVOT pre-stenting for percutaneous revaluation. For the RVOT group, CP Stent (Covered) was chosen for delivery balloon protection after rupture of the pre-dilation balloon in 7/37 patients (19%) and 30 (81%) because tear, rupture, or fracture of the conduit was expected, or further stent expansion following somatic growth was anticipated.- Sampling: n= 51 (CoA group), n=37 (RVOT group)- Mean age:<ul style="list-style-type: none">- CoA group: 19 (range from 8 to 69) years- RVOT group: 16 (range from 6 to 43) years- Sex:<ul style="list-style-type: none">- CoA group: 38M; 13F- RVOT group: 26M; 11F				P1	P2	P3		
	Report	High quality.				R1	R2	R3		
						Suitability Grade (Range 4-12)			4	
	Data Contribution	Relevant Data				Grading				

Contribution	
S&P	x
SOA	



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement


Outcomes/Endpoints	<ul style="list-style-type: none"> - Increase in diameter at coarctation (CoA group) - Decrease in peak to peak gradient (CoA group) - Number of procedures for pre-stenting and pulmonary valve delivery (RVOT Group) - increase in graft diameter (RVOT Group) - Adverse effects 	Yes 1	No 2
Follow-up	- Not specified.	Yes 1	No 2
Statistical analysis	- Two-sided $p < 0.05$ was considered significant.	Yes 1	No 2
Clinical significance	- CP Stents (Covered) can safely be applied in CHD patients. The covering allows adequate sealing of existing or expected tears, thereby increasing the safety margin with more complete dilation.	Yes 1	No 2
Data Contribution Grade (Range 4-8)		5	

Overall S&P Appraisal, Disposition and Weighting

S&P Grade (Range 9-25)	LOE (4) + Suitability (4) + Data Contribution (5) = 13	Disposition and Weighting (select)	Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21 Excluded, 22-25
---------------------------	---	------------------------------------	--

Relevant S&P Results

Safety data	<ul style="list-style-type: none"> - CoA Group: <ul style="list-style-type: none"> - No acute bleeding, aneurysm formation or life-threatening complications. - Mild procedure related-complications included groin hematoma (n = 3), transient nodal rhythm (n = 1, no wire present in left ventricle), and transient atrioventricular block with nodal escape rhythm (n = 1, while wire was present in left ventricle). - During follow-up: no stent fractures, nor stent recompression occurred, and none of the patients had limb ischemia or signs of vessel occlusion at the puncture site. - RVOT group: <ul style="list-style-type: none"> - No procedure-related complications and no extravasation. - No embolization nor fracture of CP Stent (Covered) found on annual chest X-ray follow-up.
Performance data	<ul style="list-style-type: none"> - Diameter at coarctation (CoA group): <ul style="list-style-type: none"> - Increased from 6 (0-15) to 14 (7-20) mm, $P < 0.001$. - Peak to peak gradient (CoA group): <ul style="list-style-type: none"> - Reduced from 23 (0-86) to 2 (0-25) mm Hg, $P < 0.001$. - Number of procedures for pre-stenting and pulmonary valve delivery (RVOT Group): <ul style="list-style-type: none"> - 22/37 single procedure and 15/37 in a second procedure. - Graft diameter (RVOT Group) <ul style="list-style-type: none"> - Increased from graft stenosis diameter of 13 (5-22) mm to 22 (16-26) mm at pre-revalvulation, $P < 0.001$.
Benefits/claims data	- Increase in luminal diameter in CoA patients.
Strengths	<ul style="list-style-type: none"> - CP Stent (Covered) frame is made from 90% platinum and 10% iridium 0.013" wire, welded in a zig pattern with additional gold soldering. The strut thickness is slightly larger than most other stents, but makes the stent edges relatively atraumatic. - CP Stent (Covered) was hand-crimped on a balloon-in-balloon (BIB, Numed). Hand-inflation of the balloon was performed with a 10 ml syringe on the inner balloon and 20 ml syringe on the outer balloon, automatically limiting inflation pressures to 4–6 atmospheres.
Weaknesses/	- In this retrospective study, there are no control groups with bare stents, the lack of which is inherently related to the



NuMED


Summary of Safety and Clinical Performance

SSCP – Stent Placement

		Potential bias	fact that some of these procedures would have been impossible, or significantly less safe, if bare stents were used.								
Safety & Performance (for safety only)											
Appraisal											
Level of Evidence		Study Method/Design		Question Applied			Oxford LOE 2011				
		Retrospective and observational study.		To investigate reduction in aortic wall rupture and dissection, as well as aneurysms by implanting covered stents.			1	2	3	4	5
Suitability		Relevant Data					Grading				
Device		- CP Stent (Covered) - BIB					D1	D2		D3	
Application		- CoA and re-coarctation					A1	A2		A3	
Patient		- Patients treated for CoA and re-coarctation (2 adolescents and 15 adults treated between November 2005 and January 2012). - Sampling: n=17 (11 native CoA and 6 re-coarctation) - Mean age: 35 (range 14-65) years - Sex: 4 M; 13 F					P1	P2		P3	
Report		- High quality.					R1	R2		R3	
Suitability Grade (Range 4-12)							4				
Data Contribution		Relevant Data					Grading				
Outcomes/Endpoints		- Reduction in blood pressure - Reduction in lumen diameter - Reduction of hypertensive medications at follow-up - Adverse effects					Yes 1	No 2			
Follow-up		- 2.5 years					Yes 1	No 2			
Statistical analysis		- Significance was considered as P<0.05.					Yes 1	No 2			
Clinical significance		- CP Stents (Covered) are effective in treating CoA and re-coarctation in adolescents and adults, are the treatment of choice in patients with complex anatomy, and must be available in the operating room as a rescue device when implanting a conventional stent.					Yes 1	No 2			
Data Contribution Grade (Range 4-8)							4				
Overall S&P Appraisal, Disposition and Weighting											
S&P Grade (Range 9-25)		LOE (4) + Suitability (4) + Data Contribution (4) = 12			Disposition and Weighting (select)		Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21 Excluded, 22-25				
Relevant S&P Results											
Safety data		- One death: patient died two days post-op due to massive hematemesis as a result of the combination of an extreme increase in blood pressure and an existing aneurysm. - No local complications occurred, except one hematoma that resolved spontaneously. - No patient had any complication at the iliac-femoral level that required stenting.									
Benefits/claims data		- Increased in luminal diameter - Decreased in antihypertensive medication use									

3. Alcibar et al. (2013)

Contribution	
S&P	X (safety only)
SOA	



NuMED


Summary of Safety and Clinical Performance

SSCP – Stent Placement

	Strengths	- Having observed the case of aortic rupture, and with the aim of reducing these complications in patients who have had CoA and re-coarctation since their youth, the authors decided to electively implant a NuMED (Hopkinton, New York, United States) ePTFE CP Stent (Covered). This stent is mounted on a balloon catheter and protects the vascular wall when expanded.				
	Weaknesses/ Potential bias	- Retrospective and observational study with no control group of patients receiving conventional stents. Although all patients underwent clinical follow-up, this did not include an imaging study in all cases, and so authors cannot determine with certainty the incidence of potential aneurysms.				

4. Chang et al. (2012)	Safety & Performance Appraisal										
	Level of Evidence		Study Method/Design		Question Applied		Oxford LOE 2011				
			Single arm interventional study.		To evaluate the use of CP Stent (Covered) as the primary modality in the treatment for native CoA.		1	2	3	4	5
	Suitability		Relevant Data					Grading			
	Device		- CP Stent (Covered and Bare) - BIB					D1	D2	D3	
	Application		- Native CoA					A1	A2	A3	
	Patient		- Patients with native CoA without previous treatment - Sampling: n=25 - Mean age: 22.5 (range 14-46) years - Sex: 16 M; 9 F					P1	P2	P3	
	Report		- High quality.					R1	R2	R3	
	Suitability Grade (Range 4-12)							4			
	Data Contribution		Relevant Data					Grading			
Outcomes/Endpoints		- Decrease in systolic gradient - Increase in stenotic segment diameter					Yes 1		No 2		
Follow-up		- 32 (7-72) months					Yes 1		No 2		
Statistical analysis		- P<0.05 was set as statistically significant.					Yes 1		No 2		
Clinical significance		- Implantation of CP Stent (Covered) as the primary modality is safe and effective in the treatment for native CoA in adolescents and adults. - Treatment modality of native CoA in adolescents and adults acquired excellent results, such as significant reduction in peak systolic gradient across CoA, successful relief of anatomic stenosis, and reduction of systemic hypertension. - Above all, no adverse events were encountered during the procedure or during the follow-up period of up to 72 months.					Yes 1		No 2		
Data Contribution Grade (Range 4-8)							4				
Overall S&P Appraisal, Disposition and Weighting											
S&P Grade (Range 9-25)		LOE (4) + Suitability (4) + Data Contribution (4) = 12			Disposition and Weighting (select)		Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21 Excluded, 22-25				
Relevant S&P Results											

Contribution	
S&P	x
SOA	




NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

	<table><tr><td>Safety data</td><td><ul style="list-style-type: none">- No acute complications were observed.- During a follow-up period of up to 72 months (median, 32 months and quartile range, 51 months), no adverse effects (e.g., dissection, aneurysm formation, stent migration, stent fracture) were encountered.- In the patient with the implantation of three CP stents, the aneurysm formation related to the bare CP stent was not encountered, the left subclavian artery crossed by the bare CP stent presented patent without thrombosis, and the left arm ischemia was not detected.</td></tr><tr><td>Benefits/claims data</td><td><ul style="list-style-type: none">- Reduced in peak systolic gradient.- Reduced in luminal diameter.- BIB offered precise and safe control over the stent implantation without any stent migration</td></tr><tr><td>Strengths</td><td><ul style="list-style-type: none">- Use of covered CP stents as the primary treatment modality may reduce the risk of significant complications related to stent implantation.</td></tr><tr><td>Weaknesses/ Potential bias</td><td><ul style="list-style-type: none">- Conflict of interest: not reported.</td></tr></table>	Safety data	<ul style="list-style-type: none">- No acute complications were observed.- During a follow-up period of up to 72 months (median, 32 months and quartile range, 51 months), no adverse effects (e.g., dissection, aneurysm formation, stent migration, stent fracture) were encountered.- In the patient with the implantation of three CP stents, the aneurysm formation related to the bare CP stent was not encountered, the left subclavian artery crossed by the bare CP stent presented patent without thrombosis, and the left arm ischemia was not detected.	Benefits/claims data	<ul style="list-style-type: none">- Reduced in peak systolic gradient.- Reduced in luminal diameter.- BIB offered precise and safe control over the stent implantation without any stent migration	Strengths	<ul style="list-style-type: none">- Use of covered CP stents as the primary treatment modality may reduce the risk of significant complications related to stent implantation.	Weaknesses/ Potential bias	<ul style="list-style-type: none">- Conflict of interest: not reported.		
Safety data	<ul style="list-style-type: none">- No acute complications were observed.- During a follow-up period of up to 72 months (median, 32 months and quartile range, 51 months), no adverse effects (e.g., dissection, aneurysm formation, stent migration, stent fracture) were encountered.- In the patient with the implantation of three CP stents, the aneurysm formation related to the bare CP stent was not encountered, the left subclavian artery crossed by the bare CP stent presented patent without thrombosis, and the left arm ischemia was not detected.										
Benefits/claims data	<ul style="list-style-type: none">- Reduced in peak systolic gradient.- Reduced in luminal diameter.- BIB offered precise and safe control over the stent implantation without any stent migration										
Strengths	<ul style="list-style-type: none">- Use of covered CP stents as the primary treatment modality may reduce the risk of significant complications related to stent implantation.										
Weaknesses/ Potential bias	<ul style="list-style-type: none">- Conflict of interest: not reported.										
5. Erdem et al. (2011)	Safety & Performance (for safety only)										
	Appraisal										
	Level of Evidence		Study Method/Design		Question Applied		Oxford LOE 2011				
			Single arm interventional study.		To present author's institutional experience of endovascular CP Stent implantation in children and adults with native and recurrent CoA.		1	2	3	4	5
	Suitability		Relevant Data				Grading				
	Device		<ul style="list-style-type: none">- CP Stent (16 Covered or 31 Bare) – n=47- BIB (n=29) or single balloon catheter (n=18), Z-med (not subject device)				D1		D2		D3
	Application		<ul style="list-style-type: none">- Patients with native or recurrent CoA				A1		A2		A3
	Patient		<ul style="list-style-type: none">- Patients with native CoA (Group 1); recurrent CoA and/or aneurysm developed after either surgery or balloon angioplasty (Group 2)- Sampling: n=45 (47 CP Stents, Covered or Bare)- Median age: 11 (range: 5-33) years- Sex: 34M; 11F				P1		P2		P3
	Report		<ul style="list-style-type: none">- High quality.				R1		R2		R3
			Suitability Grade (Range 4-12)				5				
Contribution											
S&P	X										
	(safety only)										
SOA											
Data Contribution		Relevant Data				Grading					
Outcomes/Endpoints		<ul style="list-style-type: none">- Decrease in invasive and echocardiographic gradients- Increase in lesion diameter- Adverse effects				Yes 1		No 2			
Follow-up		12.1±7.1 months; median 11 month (range 2-29)				Yes 1		No 2			
Statistical analysis		A p value <0.05 was considered statistically significant.				Yes 1		No 2			
Clinical significance		<ul style="list-style-type: none">- Early and short- term follow-up results indicate that stent implantation is safe and very effective in reducing coarctation gradient and increasing lesion diameter both in native and recurrent CoA.- Some serious complications do occur and hypertension remains in some patients.- Aortic disruption and stent displacement are potentially catastrophic complications of stenting				Yes 1		No 2			




NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

		but implanting a second covered stent can seal the ruptured wall and parking in a safe area or replacement of displaced stent carried by half-inflated balloon could solve the problem.						
	Data Contribution Grade (Range 4-8)				4			
	Overall S&P Appraisal, Disposition and Weighting							
	S&P Grade (Range 9-25)	LOE (4) + Suitability (5) + Data Contribution (4) = 13	Disposition and Weighting (select)	Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21 Excluded, 22-25				
	Relevant S&P Results							
	Safety data	<ul style="list-style-type: none">- No procedure related death.- Two immediate complications relating to stenting:<ul style="list-style-type: none">- One an acute wall rupture, successfully managed immediately in the same session with implantation of a second covered stent- One stent was displaced before it was completely opened. It was carried with support of partially inflated balloon and long sheath, and repositioned into the correct place.- No femoral arterial complications- No difficulty in catheter manipulation.- None of the patients required intensive care following the procedure, and all were discharged home the following day except the patient with aortic rupture and after stenting with covered stent this patient was followed two days in intensive care unit.						
	Benefits/claims data	<ul style="list-style-type: none">- Increase in luminal/lesion diameter.						
	Strengths	<ul style="list-style-type: none">- CP stent is the one of the most commonly used stent in pediatric cardiology- This stent has excellent radial strength even at larger diameters and also has brilliant visibility on fluoroscopy.						
	Weaknesses/ Potential bias	<ul style="list-style-type: none">- Some limitations have to be noted about this study:<ul style="list-style-type: none">- Firstly, there is a need a greater number of patients have undergone stent implantation and their long-term results.- Secondly, population included both children and adult.- Thirdly, this was a single-center report and patients were not compared with surgery or balloon angioplasty alone.- Fourthly, 24-hour ambulatory blood pressure monitoring before stenting was not performed in any patients.- Finally, radiologic imaging for aneurysm was done in limited number of patients after procedure.- Conflict of interest: None declared.						
	6. Butera et al. (2011)	Safety & Performance (for safety only)						
Appraisal								
Level of Evidence		Study Method/Design	Question Applied	Oxford LOE 2011				
		Prospective single arm interventional study.	To evaluate the management of aneurysms associated with CoA by covered stent deployment.	1	2	3	4	5
Contribution		Suitability		Relevant Data		Grading		
S&P	X (safety only)	Device		<ul style="list-style-type: none">- CP Stent (Covered)- BIB or Crystal balloon (not subject device)		D1	D2	D3
SOA		Application		Patients with native CoA associated with aortic wall aneurysm		A1	A2	A3
		Patient		<ul style="list-style-type: none">- Patients with CoA associated with aortic wall aneurysm- Sampling: n=11 (3 native CoA, 3 with previous surgical repair, 3 with previous balloon		P1	P2	P3



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

	angioplasty, and 2 with previous bare stent implantation)			
	- Median age: 13 (range: 6-66) years			
	- Sex: Not reported			
Report	- High quality.	R1	R2	R3
Suitability Grade (Range 4-12)		6		


Data Contribution	Relevant Data	Grading	
Outcomes/Endpoints	- Systolic pressure gradient reduction - Increase in aortic diameter - Adverse effects	Yes 1	No 2
Follow-up	- Median follow-up 50 (16-61) months	Yes 1	No 2
Statistical analysis	- P-value less than 0.05 was considered to be statistically significant	Yes 1	No 2
Clinical significance	- CP Stent (Covered) are a safe and effective treatment with low risk of complication for the treatment of CoA associated with aortic wall aneurysm. - CP Stents (Covered, e-PTFE) may be considered the treatment of choice for native CoA associated with aortic wall aneurysm.	Yes 1	No 2
Data Contribution Grade (Range 4-8)		4	

Overall S&P Appraisal, Disposition and Weighting

S&P Grade (Range 9-25)	LOE (3) + Suitability (6) + Data Contribution (4) = 13	Disposition and Weighting (select)	Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21 Excluded, 22-25
---------------------------	---	------------------------------------	--

Relevant S&P Results

Safety data	- No early complications observed.
Benefits/claims data	- Increase in luminal diameter - Reduce systolic pressure gradient - Reduce/prevent aortic wall injury (patients associated with aortic wall aneurysm)
Strengths	- Covered CP stents are manufactured with an alloy of 90% platinum and 10% iridium. Theoretically, this combination is more malleable and with good radial strength, which is enhanced by being designed in a “zig” pattern. The CP stent has rounded edges, decreasing the risk of balloon rupture or injury to the vessel wall and, in addition, the platinum component makes it more radio-opaque. Furthermore, the e-PTFE protects the stenotic and diseased segment.
Weaknesses/ Potential bias	- No conflict of interest reported.



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

7. Moltzer et al.
(2010)

Contribution	
S&P	x
SOA	

Safety & Performance

Appraisal

Level of Evidence	Study Method/Design	Question Applied	Oxford LOE 2011				
	Prospective observational study.	To evaluate the intermediate-term outcome of stent implantation for CoA in adults.	1	2	3	4	5

Suitability	Relevant Data	Grading		
Device	- CP Stent (Bare and Covered) - BIB	D1	D2	D3
Application	- Native CoA and re-coarctation	A1	A2	A3
Patient	- Patients with native CoA and re-coarctation - Sampling: n=24 - Mean age: 36 (18-60) years - Sex: 12 M; 12 F	P1	P2	P3
Report	- High quality.	R1	R2	R3
Suitability Grade (Range 4-12)			4	

Data Contribution	Relevant Data	Grading	
Outcomes/Endpoints	- Decrease in systolic gradient - Increase in minimum aortic diameter - Adverse effects	Yes 1	No 2
Follow-up	- 24 hours post intervention and 33 (8-77) months	Yes 1	No 2
Statistical analysis	- All statistical tests were two-sided and a p-value <0.05 was considered statistically significant	Yes 1	No 2
Clinical significance	- Stenting in adults results in significant blood pressure gradient decrease and increase in vessel diameter. However, serious complications do occur and hypertension remains in the majority of patients.	Yes 1	No 2
Data Contribution Grade (Range 4-8)		4	

Overall S&P Appraisal, Disposition and Weighting

S&P Grade (Range 9-25)	LOE (3) + Suitability (4) + Data Contribution (4) = 11	Disposition and Weighting (select)	Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21 Excluded, 22-25
---------------------------	---	------------------------------------	---

Relevant S&P Results

Safety data	- One death due to aorta ruptured. - Two groin hematoma post-op.
Performance data	- Systolic gradient: Decreased to < 10 mmHg in 21 patients, P<0.001 - Minimum aortic diameter: Increased from median 10 (2-17) to 16 (10-28) mm, P<0.001
Benefits/claims data	- Reduced in systolic gradient - Increased in minimum aortic diameter
Strengths	- N/A
Weaknesses/ Potential bias	- Only a small number of patients have undergone stent implantation since the authors started this procedure in 2003. This was a single-center report and patients were not compared with surgery or balloon angioplasty alone. Finally, 24-hour blood pressure monitoring before stenting was not performed in the majority of the patients. Post-stent 24-hour



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

ambulatory blood pressure monitoring is therefore difficult to translate in terms of blood pressure reduction.

8. Kische et al.
(2010)

Contribution	
S&P	x
SOA	

Safety & Performance Appraisal


Level of Evidence	Study Method/Design	Question Applied	Oxford LOE 2011				
	Technical review.	To report the technique of interventional repair in adult CoA.	1	2	3	4	5

Suitability	Relevant Data	Grading		
Device	- BIB	D1	D2	D3
Application	- CoA	A1	A2	A3
Patient	- Patients with CoA. - Sampling: Not reported. - Mean age: adult CoA patients, specific age not reported. - Sex: Not reported.	P1	P2	P3
Report	- High quality.	R1	R2	R3
Suitability Grade (Range 4-12)			5	

Data Contribution	Relevant Data	Grading	
Outcomes/Endpoints	- Stent placement (delivery of large-diameter stents). - Safety.	Yes 1	No 2
Follow-up	- Not applicable.	Yes 1	No 2
Statistical analysis	- Not applicable.	Yes 1	No 2
Clinical significance	- Not applicable.	Yes 1	No 2
Data Contribution Grade (Range 4-8)		7	

Overall S&P Appraisal, Disposition and Weighting

S&P Grade	LOE (5) + Suitability (5) +	Disposition and Weighting (select)	Accepted and Pivotal 9-12
-----------	-----------------------------	------------------------------------	---------------------------




NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

		(Range 9-25)	Data Contribution (7) = 17		Accepted but not Pivotal, 13-21 Excluded, 22-25						
		Relevant S&P Results									
		Safety data	<ul style="list-style-type: none">- BIB catheters require a larger arterial sheath for introduction, however, which needs to be upsized by 1F if a hand-crimped balloon is mounted on the balloon. Thus, although BIB catheters prevent stent flare and offer more precise control over stent placement, single-balloon catheters are still sometimes preferable in smaller patients to reduce the risk of femoral artery injury at the access site.- The use of a BIB catheter will generally prevent technical complications such as balloon rupture and stent migration.								
		Performance data	<ul style="list-style-type: none">- One of the most important technical refinements for delivery of large-diameter stents has been the NuMED Balloon-in-Balloon (BIB) catheter.- These catheters have an inner balloon and a longer outer balloon that is double the diameter of the inner balloon and are available in outer-balloon sizes of up to 24 mm.- The inner balloon of the BIB catheter is inflated, and an angiogram can be performed through the sheath or through an antegrade catheter in the proximal aorta to confirm position of the stent. With the stent in the desired position, the outer balloon is inflated to fix the stent in the lesion. Once the stent is expanded, both the outer and inner balloons are deflated as rapidly as possible.								
		Benefits/claims data	<ul style="list-style-type: none">- BIB catheters offer more precise control over stent placement								
		Strengths	<ul style="list-style-type: none">- BIB offers the important advantage of opening the stent more uniformly along its length, thereby eliminating the risk of unintended stent protrusion that has been documented by the use of single balloons.								
		Weaknesses/ Potential bias	<ul style="list-style-type: none">- No conflict of interest reported.								
9. Agnoletti et al. (2009)		Safety & Performance (for safety only)									
		Appraisal									
		Level of Evidence	Study Method/Design	Question Applied			Oxford LOE 2011				
			Two arms comparative interventional study.	To compare the CP Stent and the Palmaz stent for treatment of native and postoperative lesions of CHD patients.			1	2	3	4	5
		Suitability	Relevant Data					Grading			
		Device	<ul style="list-style-type: none">- CP Stent (Bare & Covered), crimped on BIB- Palmaz stent, crimped on BIB and simple balloons					D1	D2	D3	
		Application	<ul style="list-style-type: none">- Patients with CHD (including CoA/re-coarctation, RVOT)					A1	A2	A3	
		Patient	<ul style="list-style-type: none">- Patients with CHD (including CoA/re-coarctation, RVOT and other CHD conditions, such as transposition of the great arteries, ventricular septal defect, single ventricle, etc.)- Sampling: n= 153- 89 CP Stents (crimped on 77 BIB & 12 other balloons)- 64 Palmaz Stents (crimped on 23 BIB and 41 simple balloons)- Mean age:- CP Stents: 15.4 (SD: 9.2) years- Palmaz Stents: 11.6 (SD: 8.1) years- Sex: Not reported					P1	P2	P3	
		Report	<ul style="list-style-type: none">- High quality.					R1	R2	R3	
		Suitability Grade (Range 4-12)					6				

Contribution	
S&P	X (safety only)
SOA	



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

Data Contribution	Relevant Data	Grading	
Outcomes/Endpoints	<ul style="list-style-type: none"> - Blood pressure gradient reduction - Vessel diameter reduction - Adverse effects 	Yes 1	No 2
Follow-up	<ul style="list-style-type: none"> - Not reported. 	Yes 1	No 2
Statistical analysis	<ul style="list-style-type: none"> - A P-value less than 0.05 was considered statistically significant for stent group comparison. 	Yes 1	No 2
Clinical significance	<ul style="list-style-type: none"> - The use of the CP Stents to treat stenotic lesions of CHD is effective and relatively safe. The overall efficacy of CP Stents for the treatment of stenotic lesions is superior to that of the Palmaz stent. - CP Stents' overall safety is higher than that of the Palmaz stent; but Palmaz stents have a lower profile when inserted. 	Yes 1	No 2
Data Contribution Grade (Range 4-8)		5	

Overall S&P Appraisal, Disposition and Weighting


S&P Grade (Range 9-25)	LOE (3) + Suitability (6) + Data Contribution (5) = 14	Disposition and Weighting (select)	Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21 Excluded, 22-25
---------------------------	---	------------------------------------	--

Relevant S&P Results

Safety data	<ul style="list-style-type: none"> - Stent-related complications: <ul style="list-style-type: none"> - CP Stents (n=89 patients): 4 mild, 1 moderate, 1 severe. - Palmaz (n=64 patients): 10 mild, 2 moderate, 2 severe. - Stent migration: <ul style="list-style-type: none"> - CP Stents: 7. - Palmaz: 4. - Non stent related complications: <ul style="list-style-type: none"> - CP Stents: 1 mild, 2 moderate. - Palmaz: 1 mild, 2 moderate, 5 severe. - Urgent surgery: <ul style="list-style-type: none"> - CP Stents: 2 due to homograft rupture and stent migration. - Palmaz: 1 for aortic dissection. - Balloon related complications: Balloon burst <ul style="list-style-type: none"> - CP Stents: 0. - Palmaz: 7 (3 in BIB, 3 with simple balloons, and 1 on pre-mounted stent).
Benefits/claims data	<ul style="list-style-type: none"> - Decreased in blood pressure gradient. - Increased in vessel diameter.
Strengths	<ul style="list-style-type: none"> - Efficacy of CP Stents was similar to that of Palmaz stent for stenting of the right ventricular outflow, and higher than that of Palmaz for the stenting of aorta, but the difference was not statistically.
Weaknesses/ Potential bias	<ul style="list-style-type: none"> - Study presented retrospective results obtained in 153 consecutive patients. - CP stents were used for patients weighing more than 15 kg; and thus two populations were different concerning age, weight, indication for stenting; however, none of these differences were related to occurrence of complications. - Subgroup analyses were not performed.

10. Peters et al.

Safety & Performance



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

(2009)

Contribution	
S&P	x
SOA	

Appraisal

Level of Evidence	Study Method/Design	Question Applied	Oxford LOE 2011				
	Technical review.	To discuss the available stents and balloons in stenting in regard to their advantages and disadvantages for common applications in CHD.	1	2	3	4	5

Suitability	Relevant Data	Grading		
Device	- CP Stent and BIB	D1	D2	D3
Application	- Stenting in CoA	A1	A2	A3
Patient	- Patients with CoA - Sampling: Not reported. - Mean age: Not reported. - Sex: Not reported.	P1	P2	P3
Report	- High quality.	R1	R2	R3
Suitability Grade (Range 4-12)			6	


Data Contribution	Relevant Data	Grading	
Outcomes/Endpoints	- Design advantages or disadvantages (technical description) - Safety	Yes 1	No 2
Follow-up	- Not applicable.	Yes 1	No 2
Statistical analysis	- Not applicable.	Yes 1	No 2
Clinical significance	- Large-diameter single-balloon catheters tend to expand first at their ends and thereby evert the stent ends such that they protrude radially from the stent center. Deploying a stent in this orientation can cause injury to the vessel wall and may be a risk factor for development of aneurysm or dissection. One of the most important developments in equipment for the delivery of large-diameter stents has been the Balloon-in-Balloon (BIB; NuMED) catheter, the first balloon specifically designed for stent delivery in the CHD population.	Yes 1	No 2
Data Contribution Grade (Range 4-8)		6	

Overall S&P Appraisal, Disposition and Weighting

S&P Grade (Range 9-25)	LOE (5) + Suitability (6) + Data Contribution (6) = 17	Disposition and Weighting (select)	Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21 Excluded, 22-25
---------------------------	---	------------------------------------	--

Relevant S&P Results

Safety data	<ul style="list-style-type: none"> - BIB: - Stent foreshortening can be much higher if stents are expanded with a single large balloon directly to 15 or 18 mm. The reason for this is that the ends of the stent are compressed toward each other due to the typical “dumb-belling” of the balloon at the end of inflation while the center of the stent is expanding to its full diameter. This results in significant shrinkage of the overall length of the stent. Thus, if final length is critical, the delivery requires sequential balloon dilatation with increasing diameters or, even better, the use of a Balloon-in-Balloon (BIB™ catheter; NuMED). - Balloon rupture with inadequate stent expansion may be prevented by avoiding kinking of the balloon/stent assembly by the use of newer stents with softer ends and by the use of BIB systems.
-------------	---



NuMED


Summary of Safety and Clinical Performance

SSCP – Stent Placement

	Performance data	<div>- BIB:<ul style="list-style-type: none">- These catheters have an inner balloon and a longer outer balloon that is double the diameter of the inner balloon. The BIB catheters offer the important advantage of opening the stent more uniformly along its length but require a larger arterial sheath for introduction.- While BIB catheters prevent stent flare and offer more precise control over stent placement, single-balloon catheters are still sometimes preferable in smaller patients to reduce risk of injury to the femoral artery at the access site.</div>									
	Benefits/claims data	<div>- BIB offers more precise control over stent placement</div>									
	Strengths	<div>- CP Stent:<ul style="list-style-type: none">- These stents have excellent visibility on fluoroscopy and maintain excellent radial strength even at larger diameters.</div>									
	Weaknesses/ Potential bias	<div>- No conflict of interest reported.</div>									

11. Tzifa et al. (2006)	Safety & Performance Appraisal											
	Level of Evidence		Study Method/Design		Question Applied			Oxford LOE 2011				
			Single arm interventional study.		To evaluate the use of Covered CP Stents in treatment of CoA.			1	2	3	4	5
	Suitability		Relevant Data					Grading				
	Device		<div>- CP Stent (Covered) - BIB</div>					D1	D2		D3	
	Application		<div>- CoA</div>					A1	A2		A3	
	Patient		<div>- Patients with CoA (fully grown patients) - Sampling: n=30 - Mean age: 28±17.5 (range 8 to 65) years - Sex: not reported</div>					P1	P2		P3	
	Report		<div>- High quality.</div>					R1	R2		R3	
	Suitability Grade (Range 4-12)							5				
	Data Contribution		Relevant Data					Grading				
	Outcomes/Endpoints		<div>- Reduction in blood pressure gradient - Reduction in coarctation diameter</div>					Yes 1		No 2		
	Follow-up		<div>- 11 months</div>					Yes 1		No 2		
	Statistical analysis		<div>- Statistical significance was defined as P<0.05.</div>					Yes 1		No 2		
	Clinical significance		<div>- CP Stents (Covered) may be used as the therapy of choice in patients with complications after CoA repairs, whereas they provide a safe alternative to conventional stenting in patients with severe and complex CoA lesions or advanced age.</div>					Yes 1		No 2		
	Data Contribution Grade (Range 4-8)							4				
Overall S&P Appraisal, Disposition and Weighting												
S&P Grade (Range 9-25)		LOE (4) + Suitability (5) + Data Contribution (4) = 13			Disposition and Weighting (select)		Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21					

Contribution	
S&P	x
SOA	



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

Excluded, 22-25

Relevant S&P Results

Safety data	- Two stent fractures in the “old” design of the stent, no fractures in the “new” stent design Note: Since May 2002, the CP Stents (Covered) have been produced with reinforced golden soldering joints as the “new” stent design
Performance data	- Blood pressure gradient: From 36 + 20 mmHg to 4 + 4 mmHg, P<0.0001 - Diameter at coarctation: From 6.4 +3.8 mm to 17.1 + 3.1 mm, P<0.0001
Benefits/claims data	- Reduction in blood pressure gradient - Reduction in coarctation diameter - BIB allows readjustment of position after inflation of the inner balloon.
Strengths	- Covered stents were chosen: 1) as a rescue treatment in patients with CoA aneurysms or previous stent-related complications; and 2) in patients at risk of complications because of complex CoA anatomy or advanced age (defined as >65 years) - Covered CP stents are made of a framework of platinum iridium wire welded in a zig pattern. The addition of a gold soldering to each weld spot fills any voids caused by the welding and transfers the stresses to a larger area of the stent. The gold also serves to encapsulate the welded area, once again adding to the total strength of the weld. The stent is then fitted with a covering of ePTFE to achieve a solid tubular structure that retains fluid. The ePTFE covering is initially approximately 7 mm in diameter and will stretch over the range of diameters of expansion (usually from 12 to 24 mm diameter), and will always be taut over the stent when expanded. When the covering is mounted, it is folded over the crimped stent and expands uniformly when the balloon is inflated. - The BIB allows for readjustment of position after inflation of the inner balloon.
Weaknesses/ Potential bias	- Not reported.


Safety & Performance Appraisal

Level of Evidence	Study Method/Design	Question Applied	Oxford LOE 2011				
	Comparative, two single arm interventional study (CP Stent/BIB versus Palmaz stent/single balloon).	To demonstrate effectiveness of CP Stent, in combination with BIB, for treating aortic coarctation in comparison with the Palmaz stent.	1	2	3	4	5

12. Cheatham et al. (2001a)

Contribution	
S&P	x
SOA	

Suitability	Relevant Data	Grading		
Device	- CP Stent and BIB. Note: D1 for subject device (CP Stent and BIB). - Palmaz stent and single balloon.	D1	D2	D3
Application	- CoA and re-coarctation	A1	A2	A3
Patient	- Patients with CoA and re-coarctation - Sampling: n=46 (21 Palmaz Stent, 25 CP Stent) - Mean age: - Palmaz Stent: 12 (range: 4.5 to 16) years old - CP Stent: 24.1 (range: 10.5 to 60) years old - Sex: M:F - Palmaz Stent: 15M; 6F - CP Stent: 15M; 10F	P1	P2	P3
Report	- High quality.	R1	R2	R3



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement


Suitability Grade (Range 4-12)			4
Data Contribution	Relevant Data	Grading	
Outcomes/Endpoints	<ul style="list-style-type: none"> - Decrease in peak systolic gradient - Safety 	Yes 1	No 2
Follow-up	<ul style="list-style-type: none"> - Limited follow-up because of the relatively short-elapsd time and multiple, out of country institutions involved 	Yes 1	No 2
Statistical analysis	<ul style="list-style-type: none"> - Statistically significant achieved when $P < 0.05$. 	Yes 1	No 2
Clinical significance	<ul style="list-style-type: none"> - Both the Palmaz and the NuMED CP Stents offer an effective, nonsurgical treatment for both native and recurrent CoA, regardless of site or severity of obstruction. - Native CoA tends to be more severe, with tighter stenosis and higher gradients compared to recurrent coarctation. Aneurysm development may occur in these patients after Palmaz stent implantation. Therefore, graduated serial stent dilation and/or covered stent implantation should be considered in high-risk patients. 	Yes 1	No 2
Data Contribution Grade (Range 4-8)			5

Overall S&P Appraisal, Disposition and Weighting

S&P Grade (Range 9-25)	LOE (3) + Suitability (4) + Data Contribution (5) = 12	Disposition and Weighting (select)	Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21 Excluded, 22-25
---------------------------	---	------------------------------------	--

Relevant S&P Results

Safety data	<ul style="list-style-type: none"> - Intraoperative complications: <ul style="list-style-type: none"> - Two cases of stent embolization in Palmaz Stent group versus one case of left hemothorax in CP Stent (Covered) group - Late complications: <ul style="list-style-type: none"> - Aortic aneurysm in Palmaz Stent Group (n=3) - 9/21 (43%) continued to require antihypertensive medication - Adverse effects: <ul style="list-style-type: none"> - Flaring of the Palmaz Stent as a result of the single balloon that first expands at the ends, not from the middle. - Premature BIB rupture reported using Palmaz stent, believed to be due to inadvertent puncture of outer balloon by the sharp-edged Palmaz stent during hand crimping.
Performance data	<ul style="list-style-type: none"> - Peak systolic gradient <ul style="list-style-type: none"> - Palmaz Stent <ul style="list-style-type: none"> - Native coarctation: 46.8 to 1.5 mmHg ($P < 0.05$) - Recurrent coarctation: 35 to 1.2 mmHg ($P < 0.05$) - CP Stent <ul style="list-style-type: none"> - Native coarctation of the aorta: 53 mmHg to 2 mmHg ($P < 0.001$) - Recurrent coarctation: 41 mmHg to 1.2 mmHg ($P < 0.001$)
Benefits/claims data	<ul style="list-style-type: none"> - BIB significantly improves stent delivery and final deployment of any stent. The inner balloon is always inflated first, partially expanding the stent without flaring and allows repositioning of the stent before final deployment, when the outer balloon is inflated.
Strengths	<ul style="list-style-type: none"> - The zig design of the NuMED CP Stent improved strength and flexibility while minimizing stent shortening and vessel/balloon trauma.



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

		<ul style="list-style-type: none">- The tempered platinum/iridium wire and zig design of the NuMED CP Stent improved strength, flexibility, and radiopacity while minimizing stent shortening and vessel/balloon trauma.- The NuMED CP stent offers a wider range of expanded diameters and lengths than the Palmaz stent, which is necessary to ensure an adequate adult vessel diameter in growing children and young adults and reduces the need for multiple stents in long segment obstruction.- The BIB catheter is selected so that the inner balloon is always shorter than the stent while the outer balloon is slightly longer. The BIB delivery catheter significantly improves stent delivery and final deployment of any stent while minimizing stent migration and flaring. This in turn decreases the incidence of ventricular tachycardia in interventional pediatric cardiologists.
	Weaknesses/ Potential bias	N/A

13. Cheatham et al. (2001b)	Safety & Performance Appraisal									
	Level of Evidence		Study Method/Design		Question Applied			Oxford LOE 2011		
			Report results collected 45 patients underwent CP stent implantation from August 1998 through August 1999.		To report the CP Stent and BIB development, including results from clinical study conducted from August 1998-August 1999.			1	2	3
	Suitability		Relevant Data					Grading		
	Device		<ul style="list-style-type: none">- CP Stent- BIB					D1D2D3		
	Application		<ul style="list-style-type: none">- CoA and other conditions.					A1A2A3		
	Patient		<ul style="list-style-type: none">- Patients with CoA and other conditions:<ul style="list-style-type: none">- CoA (n=25; 17 native CoA and 8 Re-coarctation),- Right pulmonary artery (RPA) stenosis (n=5),- Isolated left pulmonary artery (LPA) stenosis (n=2),- Bilateral branch stenosis (n=6),- Recurrent right ventricle to pulmonary artery (RV-PA) homograft stenosis (n=4),- Blalock-Taussig shunt stenosis (n=1),- Multiple sites of left-to-right shunt inside a lateral tunnel Fontan repair (n=1),- Obstructed superior vena cava (SVC) baffle limb after Mustard repair for transposition of the great arteries (n=1)- Sampling: n=45 patients (CP Stent, n=57)- Mean age: 19 (range 1.8-60) years- Sex: 25 M; 20 F					P1P2P3		
	Report		<ul style="list-style-type: none">- High quality.					R1R2R3		
	Suitability Grade (Range 4-12)								6	
	Data Contribution		Relevant Data					Grading		
	Outcomes/Endpoints		<ul style="list-style-type: none">- Peak systolic gradient reduction- Procedural complications					Yes 1No 2		
	Follow-up		<ul style="list-style-type: none">- Minimal follow-up due to short study period and large number of institutions involved					Yes 1No 2		
	Statistical analysis		<ul style="list-style-type: none">- P-values reported.					Yes 1No 2		

Contribution	
S&P	x
SOA	



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement


Clinical significance	<ul style="list-style-type: none"> - NuMED CP stent (placed by BIB stent placement catheter) offers an effective, non-surgical treatment for a wide variety of vascular obstructions associated with congenital heart disease. - The NuMED BIB catheter is an innovative concept that has significantly improved operator control during intravascular stent delivery. It minimizes stent flaring, migration, and shortening, while effectively eliminating catheter movement during deployment. It also allows the partially expanded stent to be repositioned before final expansion, which is a significant benefit to the interventionalist to maintain control and precisely position the stent. We currently use the BIB catheter for all stents. Although more data and longer follow-up are required, the NuMED CP stent and BIB delivery catheter offer great promise in the future treatment of children and adults with congenital heart disease. 	Yes 1	No 2
Data Contribution Grade (Range 4-8)		5	

Overall S&P Appraisal, Disposition and Weighting

S&P Grade (Range 9-25)	LOE (4) + Suitability (6) + Data Contribution (5) = 15	Disposition and Weighting (select)	Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21 Excluded, 22-25
---------------------------	---	------------------------------------	--

Relevant S&P Results

Safety data	<ul style="list-style-type: none"> - Two procedural complications, both considered avoidable and there were no further sequelae - One with severe native CoA requiring a covered stent, there was transient left hemothorax. - One traumatic stent fracture during attempted entry of the long, covered CP stent in the Fontan patient using a modified 'front-load' technique. The stent was inadvertently pushed out of the delivery sheath in the groin with the first row of 10 zigs being traumatized and fractured. - Follow-up: <ul style="list-style-type: none"> - Stent fatigue fracture and fragment embolization in two patients - Two patients with severe native CoA and stenoses <2 mm had immediate residual gradients of 20 and 25 mmHg secondary to limited stent expansion to avoid excessive vessel trauma and possible aneurysm formation with planned stent re-dilation later. - One patient had a 30 mmHg residual aortic gradient 10 months post implant secondary to an intimal flap that was successfully treated with a second CP stent
Performance data	<ul style="list-style-type: none"> - Peak systolic gradient was reduced in 17 patients with native CoA from 56.2 mmHg to 4.6 mmHg, while in the 8 patients with recurrent aortic obstruction, the gradient was reduced from 41.8 mmHg to 0.9 mmHg, both statistically significant at P<0.001 using paired t-tests - Isolated RPA and LPA stenoses were also effectively treated with peak systolic gradient reductions from 54.6 to 5 mmHg and 52.5 to 6.5 mmHg respectively P<0.001 - In the six children with combined RPA and LPA stenoses, 'kissing stents' reduced the peak systolic gradients from 43.5 and 45 mmHg to 6.8 and 6.0 mmHg, respectively P<0.01. - The four patients with recurrent RV-PA homograft obstruction also had effective relief of their gradients from 55 to 14.3 mmHg P<0.01. - After stenting the stenotic right Blalock-Taussig shunt in the young man with complex cyanotic congenital heart disease, O2 saturations increased from 78 to 88%. The implantation of the long covered CP stent was also clinically effective in treating the young man with multiple leaks in the lateral tunnel by improving resting O2 saturations from 80 to 96%. - Finally, the 9 mmHg mean gradient across the obstructed SVC baffle after Mustard's repair was completely eliminated.



NuMED


Summary of Safety and Clinical Performance

SSCP – Stent Placement

	Benefits/claims data	<ul style="list-style-type: none">- The BIB effectively eliminating catheter movement during deployment. It also allows the partially expanded stent to be repositioned before final expansion, which is a significant benefit to the interventionalist to maintain control and precisely position the stent.									
	Strengths	<ul style="list-style-type: none">- CP Stent versus Palmaz Stent<ul style="list-style-type: none">- The advantages of the NuMED CP stent compared to the Palmaz stent are as follows: (1) superior radiopacity secondary to the platinum composition; (2) superior ‘compression’ or radial hoop strength secondary to the tempered wire and zig design; (3) less rigidity because of the malleability of the tempered platinum-iridium wire; (4) less potential trauma to the delivery balloon and target vessel secondary to the rounded edges of the zig pattern; (5) wider range of expanded diameters from 8 to 24 mm (10 zig can be expanded to 30 mm) while maintaining ≤20% stent shortening; (6) superior selection of stent lengths to meet the demands of a wide range of target lesions; and (7) maximal stent shortening of <20% will minimize chances of missing the target site or need of multiple serial stents.- BIB versus single balloon catheter:<ul style="list-style-type: none">- Careful observation of how a single balloon catheter may actually create significant problems during any stent deployment, but is exaggerated in the aorta. Conventionally, the balloon is chosen to be longer than the stent to avoid stent migration during delivery or deployment. Unfortunately, the proximal and distal ends of the balloon catheter expand first and well before the stent. This leads to flaring of the edges of the stent, which in the Palmaz stent’s case is dangerous because of the sharp leading and trailing edges approaching the vessel wall and balloon leading to trauma of both. In addition, the partially expanded balloon acts as a floatation catheter, allowing catheter and stent movement prematurely during deployment. Finally, there is no ability to reposition the stent during deployment.- In November 1997, the NuMED Balloon In Balloon (BIB) catheter was designed with an inner Tyshak balloon and an outer Z-Med balloon. The inner balloon is very low profiled and expands to half the outer balloon diameter, while the length is 1 cm shorter than the outer balloon. The inner balloon is always inflated first using a twisting action of the locked endoflator that expands the stent to 0.5 of the target vessel diameter without flaring of ends of the stent, since the balloon is shorter than the stent. Because the stent is still in contact with the unexpanded outer balloon material, the entire stent-balloon delivery catheter system can be repositioned before final deployment by expanding the outer balloon.									
	Weaknesses/ Potential bias	A Tower, DJ Villnave and R Normile (NuMED) provided technical support for this publication.									

14. Meadows et al. (2015)	Safety & Performance This publication presents the results from the COAST trial for CP Stent (Bare and Covered) to treat native and recurrent CoA in selected children, adolescents and adult (NCT00552812). Please refer information presented in Table G-1 for safety and performance of the subject devices, Study no. 1.									
	Contribution									
	S&P	x								
SOA										

15. Bairam et al. (2021)	Safety & Performance Appraisal											
	Level of Evidence	Study Method/Design	Question Applied				Oxford LOE 2011					
		Single center prospective study.	The aim of this study was to evaluate the immediate and intermediate results of CP stenting for native CoA.				1	2	3	4	5	
Contribution		Suitability				Relevant Data				Grading		
S&P	x	Device		- Covered CP stent was hand crimped on a balloon dilation catheter (either Z-med balloon				D1	D2	D3		
SOA	x											



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

	dilatation catheter or BIB balloon according to availability in the stock)			
	- A total of 39 stents implanted; the CP stent used in 30 patients, and uncovered stent was used in one patient to avoid occlusion of the left subclavian artery			
Application	- Native CoA	A1	A2	A3
Patient	- 32 consecutive adult patients with native CoA who underwent aortic stent replacement - Referred with clinical evidence for CoA (hypertension, arm-leg blood pressure difference ≥ 20 mmHg), for whom the diagnosis of CoA was confirmed by both transthoracic echocardiography and by CT angiography - Mean age 30.83 ± 11.179 years (16-56) - 15 (46.9%) male and 17 (53.1%) female	P1	P2	P3
Report	- Publication reports the CP stent was implanted in 30 patients and uncovered stent in one patient but does not report the total of each stent design implanted.	R1	R2	R3
Suitability Grade (Range 4-12)		5		


Data Contribution	Relevant Data	Grading	
Outcomes/Endpoints	- Short-term pre/post-implant hemodynamics and angiographic data were reported. Clinical, echocardiographic and CT exam for restenosis, aneurysm formation, stent migration, changes in blood pressure, and the use of antihypertensive drugs were recorded during follow-up.	Yes 1	No 2
Follow-up	- Ranged from 6-32 months (mean 14.9 months) included clinical, echocardiographic and CT exam assessing for restenosis, aneurysm formation, stent migration, blood pressure control and the intensity of antihypertensive medications	Yes 1	No 2
Statistical analysis	- Mean, standard deviation and P-values reported	Yes 1	No 2
Clinical significance	- The magnitude of the treatment effect observed was clinically significant.	Yes 1	No 2
Data Contribution Grade (Range 4-8)		4	

Overall S&P Appraisal, Disposition and Weighting

S&P Grade (Range 9-25)	LOE (3) + Suitability (5) + Data Contribution (4) = 12	Disposition and Weighting (select)	Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21 Excluded, 22-25
---------------------------	---	------------------------------------	--

Relevant S&P Results

Safety data	- There was no major complications, with no deaths. - There was no recorded complication at the implantation site except for one case of ruptured balloon used to deliver a second CP stent, the balloon retrieved to the iliac artery then successfully removed surgically. - No significant complications were seen during procedure and at six months follow up.
Performance data	- Procedure success rate was 93.4% (n=31). The technique was considered effective if the invasive grade was decrease to < 20 mmHg and increased the angiographic diameter $> 50\%$. - Peak gradient across the coarctation fell from 60.0 ± 21.960 to 10.0 ± 19.821 mm Hg post procedure ($P = 0.0001$). - Systolic blood pressures fell from 164.6 ± 25.889 mm Hg to 138.1 ± 17.006 mm Hg immediately after stenting and 134.3 ± 12 mm Hg at six months.
Benefits/claims data	- "...BIB balloon is an good instrument to prevent stent migration and rupture of balloon during the technique." - Treating native CoA with stent implantation is a less invasive, safe and highly successful technique with excellent angiographic and clinical response that is sustained to the midterm follow up of this study.



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

Strengths	<ul style="list-style-type: none"> - Prospective design - Consecutive patients
Weaknesses/ Potential bias	<ul style="list-style-type: none"> - Single center design - Follow up ranged from 6-32 months (mean 14.9 Months) although study intended as immediate and intermediate results

State of the Art

Appraisal

Medical condition	Alternatives	Risk/benefit	Side-effects	Equivalence	Surrogate endpoints
Yes 1 No 2	Yes 1 No 2	Yes 1 No 2	Yes 1 No 2	Yes 1 No 2	Yes 1 No 2

Overall SOA Appraisal and Disposition

SOA Grade (Range 6-12)	10	Disposition (select)	Accepted, < 12 Excluded, 12
---------------------------	----	----------------------	--------------------------------

Relevant SOA Results

SOA data	Risk <ul style="list-style-type: none"> - There was no recorded complication at the implantation site except for one case of ruptured balloon used to deliver a second CP stent, the balloon retrieved to the iliac artery then successfully removed surgically.
Comments	- No author-identified limitations. Funding and conflicts of interest not addressed.

Safety & Performance


Appraisal

Level of Evidence	Study Method/Design	Question Applied	Oxford LOE 2011				
	Single center retrospective study.	The aim of this study was to analyze long-term results after CoA treatment with bare and covered CP stents in our institution and to derive recommendations for the differential use of these stent types.	1	2	3	4	5

Suitability	Relevant Data	Grading		
Device	<ul style="list-style-type: none"> - 212 patients received treatment with bare (n=71) and covered (n=141) CP stents between September 1999 and July 2021 - Stents were mounted on BIB catheters, not specified as pre-mounted 	D1	D2	D3
Application	- Native CoA (n=110/212, 51.9%) and recoarctation after primary surgical or interventional treatment (n=102/212, 48.1%)	A1	A2	A3
Patient	<ul style="list-style-type: none"> - Median patient age was 18.8 years (IQR 11.9; 35.8) - Median patient weight 61.3 kg (IQR 43.3; 74.7) - 146/212 (68.9%) male 	P1	P2	P3
Report	- High quality.	R1	R2	R3
Suitability Grade (Range 4-12)		4		

16. Schleiger et al.
(2023)

Contribution	
S&P	x
SOA	x



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement


Data Contribution	Relevant Data	Grading	
Outcomes/Endpoints	- Procedural success, survival rate, freedom from re-intervention, peri-procedural and long-term complications were reported.	Yes 1	No 2
Follow-up	- Median follow-up of 7.3 years (IQR: 4.3-12.6)	Yes 1	No 2
Statistical analysis	- Patient characteristics expressed as median and IQR. Survival and freedom from re-intervention were assessed using Kaplan-Meier survival analysis. Survival and reintervention rates between groups were compared using the log rank test. Differences between groups were analyzed using the χ^2 test for categorical variables and Wilcoxon rank sum test for continuous variables. Potential risk factors for re-intervention were evaluated with univariate logistic and Cox regression analysis. Time-independent variables were included in a multivariable model using HR. A $p < 0.05$ was considered statistically significant.	Yes 1	No 2
Clinical significance	- The magnitude of the treatment effect observed was clinically significant.	Yes 1	No 2
Data Contribution Grade (Range 4-8)		5	

Overall S&P Appraisal, Disposition and Weighting

S&P Grade (Range 9-25)	LOE (3) + Suitability (4) + Data Contribution (4) = 11	Disposition and Weighting (select)	Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21 Excluded, 22-25
---------------------------	---	------------------------------------	---

Relevant S&P Results

Safety data	<ul style="list-style-type: none"> - Survival rate: Survival rate was 98.1% after five, and 95.6% after 10 and 15 years, respectively, and did not differ between patients who received bare or covered CP stents (Log Rank $p = 0.263$). In-hospital mortality occurred in 1/212 patients (0.5%) and late mortality in 8/158 patients (5.1%). Late mortality was not attributable to previous CoA treatment...There was no difference in late mortality according to stent type ($p = 0.261$). - Complications rate: <ul style="list-style-type: none"> o Peri-procedural complications – Entire cohort ($n=212$) <ul style="list-style-type: none"> ▪ Injury/thrombosis of vascular access vessel: 9/212 (4.2%) ▪ Bleeding of vascular access vessel: 1/212 (0.5%) ▪ Stent dislocation: 2/212 (0.9%) <ul style="list-style-type: none"> • In two patients with an unsuccessful procedure the stent migrated into the descending aorta immediately after placement. (interpreted as stent migration) ▪ Aortic dissection/aortic wall rupture: 3/212 (1.4%) o Long-term complications – Entire cohort <ul style="list-style-type: none"> ▪ Aneurysm formation: 14/133 (10.5%) ▪ Stent fracture: 19/108 (17.6%) - In our study, no association between anatomic and hemodynamic characteristics (e.g., minimal CoA diameter, maximal peak-to-peak systolic gradient) or technical procedural details (BiB® catheter diameter, achieved nominal pressure, postdilatation) and the occurrence of aortic wall injuries was detected.
Performance data	<ul style="list-style-type: none"> - Procedural success was achieved in 187/212 (88.2%) patients. - After stent implantation a significant reduction of systolic blood pressure was achieved from a preinterventional median pressure of 145 mmHg (IQR 134; 157) to a postinterventional median pressure of 123 mmHg (IQR 112; 135) ($p < 0.001$). - In 25 patients the interventional procedure was not considered successful: In 22 of these patients a reduction of the peak systolic pressure <10 mmHg was not achieved.



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

	<ul style="list-style-type: none">○ Fifteen of these patients were additionally diagnosed with a hypoplastic aortic arch with a remaining systolic ascending to descending aortic pressure difference >10 mmHg after successful implantation of the CP stent in the CoA region.○ In two patients with an unsuccessful procedure the stent migrated into the descending aorta immediately after placement. Both patients underwent subsequent surgical repair after fixation of the stent in the descending aorta by balloon dilatation.○ Another patient developed an aortic wall rupture immediately after stent implantation and died during extracorporeal cardiopulmonary resuscitation. <ul style="list-style-type: none">- Re-intervention rate: Planned re-interventions were performed in 33/158 patients (20.9%)...In 44/158 patients (27.8%), unplanned re-interventions were performed to treat re-stenosis or aortic wall injuries... The probability of freedom from re-intervention in the entire cohort was 81.0% after five, 64.0% after 10 and 62.0% after 15 years, respectively. Re-intervention rate did not differ between patients who received endovascular CoA treatment with bare or covered CP stents ($p = 0.50$)... Multivariable risk factor analysis revealed previous CoA surgery (HR: 2.0, 95% confidence interval [CI]: 1.1–3.9, $p = 0.029$), postdilatation (HR: 2.9, 95% CI: 1.1–6.3, $p = 0.028$) and age at intervention (HR: 0.96, 95% CI: 0.94–0.99, $p = 0.002$) as independent risk factors for re-intervention.- Antihypertensive medications: Residual arterial hypertension was present in 53/158 patients (33.5%); 23 of these patients were adult and 30 patients pediatric ($p = 0.173$). Before endovascular treatment with bare or covered CP stent 114 patients received no medical antihypertensive therapy, 44 patients received monotherapy, 30 patients dual therapy, and 24 patients triple or quadruple therapy. In 76/158 patients (48.1%) the number of antihypertensive medications did not change during follow-up, whereas in 57/158 patients (36.1%) the number of antihypertensive medication was increased and in 25/158 patients (15.8%) decreased.- Median diameter of BIB catheter (entire cohort, $n=212$): 16.0mm [IQR: 14.0-20.0]- Full expansion of BIB catheter (entire cohort, $n=212$): 171/212 (80.7%)
Benefits/claims data	<ul style="list-style-type: none">- In conclusion, our study documents excellent long-term results after CoA treatment with bare and covered CP stents. In our cohort, mortality, re-intervention and complication rate did not significantly differ between both stent types.
Strengths	<ul style="list-style-type: none">- Long-term follow-up (median of 7.3 years, IQR: 4.3-12.6)
Weaknesses/ Potential bias	<ul style="list-style-type: none">- Single center design- Retrospective design

State of the Art

Appraisal


Medical condition		Alternatives		Risk/benefit		Side-effects		Equivalence		Surrogate endpoints	
Yes 1	No 2	Yes 1	No 2	Yes 1	No 2	Yes 1	No 2	Yes 1	No 2	Yes 1	No 2

Overall SOA Appraisal and Disposition

SOA Grade (Range 6-12)	10	Disposition (select)	Accepted, < 12 Excluded, 12
---------------------------	----	----------------------	--------------------------------

Relevant SOA Results

SOA data	<p>Complications rate:</p> <ul style="list-style-type: none">- Peri-procedural complications – Entire cohort ($n=212$)<ul style="list-style-type: none">○ Injury/thrombosis of vascular access vessel: 9/212 (4.2%)○ Bleeding of vascular access vessel: 1/212 (0.5%)○ Stent dislocation: 2/212 (0.9%)
----------	---



NuMED


Summary of Safety and Clinical Performance

SSCP – Stent Placement

		<ul style="list-style-type: none">○ Aortic dissection/aortic wall rupture: 3/212 (1.4%)- Long-term complications – Entire cohort<ul style="list-style-type: none">○ Aneurysm formation: 14/133 (10.5%)○ Stent fracture: 19/108 (17.6%)				
	Comments	<ul style="list-style-type: none">- Author-identified limitations included unequal distribution between subgroups with more patients receiving Covered CP stents due to institutional preference, difference in follow-up duration between sub-groups due to differences in time points of availability and product approval, 54 patients lost to follow-up could not be considered in long-term analysis, the low event rate of major complications (e.g., aortic dissection, stent fracture, aneurysm formation) may limit statistical comparison between groups, incidence of long-term aortic wall complications may be underestimated as re-catheterization or cross-sectional imaging was not available for all patients, non-invasive blood measurement during exercise or 24-hours blood pressure measurements were not available to identify unmasked arterial hypertension and the number of antihypertensive medications may be affected by cofounders.- Authors declare no conflict of interest. Open access funding enabled and organized by Projekt DEAL.				

17. Vargas-Acevedo et al. (2024)	Safety & Performance Appraisal								
	Level of Evidence		Study Method/Design		Question Applied		Oxford LOE 2011		
			Single-center retrospective review		Treatment Benefits, Treatment Harms		1	2 3 4 5	
	Suitability		Relevant Data			Grading			
	Device		<ul style="list-style-type: none">- Covered CP Stent- BIB			D1	D2	D3	
	Application		<ul style="list-style-type: none">- CoA and variants of arch obstruction- 76% (n=19) mechanically crimped, 24% (n=6) manually crimped			A1	A2	A3	
	Patient		<ul style="list-style-type: none">- Sampling: n=25 eligible events- Median age: 18 years (IQR: 18-28)- Sex: 10 M; 15 F- Diagnosis: Isolated CoA (80%), CoA with hypoplastic aortic arch (12%), interrupted aortic arch (8%)			P1	P2	P3	
	Report		<ul style="list-style-type: none">- High quality			R1	R2	R3	
	Suitability Grade (Range 4-12)						6		
	Data Contribution		Relevant Data			Grading			
	Outcomes/Endpoints		<ul style="list-style-type: none">- Median sheath size			Yes 1	No 2		
	Follow-up		<ul style="list-style-type: none">- None reported			Yes 1	No 2		
	Statistical analysis		<ul style="list-style-type: none">- P-values reported			Yes 1	No 2		
	Clinical significance		<ul style="list-style-type: none">- Device delivery through smaller sheaths could permit interventions to be performed in smaller patients while reducing the risk of significant vascular access related injury without compromising procedure or stent performance.			Yes 1	No 2		
	Data Contribution Grade (Range 4-8)						5		
	Overall S&P Appraisal, Disposition and Weighting								
	S&P Grade (Range 9-25)		LOE (3) + Suitability (6) + Data Contribution (5) = 14			Disposition and Weighting (select)		Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21	

Contribution	
S&P	X (safety only)
SOA	



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

Excluded, 22-25

Relevant S&P Results

Safety data	- There was no record of discarded stents, no balloon ruptures, nor need to upsize the sheath in the mechanically crimped cohort.
Performance data	- N/A
Benefits/claims data	- N/A
Strengths	- The median sheath size for mechanically crimped stents was - 2 Fr compared with a median of 0 Fr for manually crimped stents ($p = 0.007$).
Weaknesses/ Potential bias	<ul style="list-style-type: none"> - Author-identified limitations include lack of routine record of crimping method, details of difficulties with stent behavior during manual crimping and advancement through the delivery sheath in many procedural notes, limited population size for retrospective analysis, non-randomized single center experience which can limit the power of the results. - No funding was used in this study. - Dr. Morgan is a consultant for NUMED Inc. The remaining authors declared no conflict of interest.

Safety & Performance Appraisal


Level of Evidence	Study Method/Design	Question Applied	Oxford LOE 2011				
	Open-label, parallel-group, blinded endpoint randomized pilot clinical trial.	To compare the safety and efficacy of the balloon-expandable stent (BES) and the self-expandable stent (SES) in the endovascular treatment of CoA.	1	2	3	4	5

Suitability	Relevant Data	Grading		
Device	<ul style="list-style-type: none"> - Uncovered CP BES (n=46 patients) - Uncovered nitinol SES (n=46 patients) 	D1	D2	D3
Application	- Adult patients with de novo native CoA	A1	A2	A3
Patient	<ul style="list-style-type: none"> - 92 eligible patients (32 women, 34.8%) with a median age of 30 years (IQR: 20-36 years) were randomized equally into the two groups - Median age <ul style="list-style-type: none"> o BES: 29.9 years (IQR: 19.5-37.0 years) o SES: 28.6 years (IQR: 21.0-33.5 years) - Female <ul style="list-style-type: none"> o BES: 14 (30.4%) o SES: 18 (39.1%) 	P1	P2	P3
Report	- Report suitable for review.	R1	R2	R3
Suitability Grade (Range 4-12)		4		

Data Contribution	Relevant Data	Grading	
Outcomes/Endpoints	<ul style="list-style-type: none"> - Primary outcome was a composite of periprocedural and vascular access complications. - Secondary outcomes were composed of the incidence of aortic recoarctation, thoracic aortic aneurysm/pseudoaneurysm, and residual hypertension at the 12-month follow-up. 	Yes 1	No 2
Follow-up	- 12-month follow-up period at intervals of one, three, six and 12 months	Yes 1	No 2
Statistical analysis	- Data were expressed as mean \pm SD or median (IQR) for interval variables and counts (%) for	Yes 1	No 2

18. Sadeghipour et al. (2022)

Contribution	
S&P	x
SOA	x



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

	<p>categoric variables. Categoric variables were compared using the chi-square test or the Fisher exact test. Continuous variables were compared between the two groups with the aid of the Student's t-test (or its nonparametric equivalent, the Mann-Whitney U test). Analyses of the study outcomes were based on the binary logistic regression and the odds ratio (OR) as the effect size. The cumulative incidence rate of the primary composite endpoint with its respective 95% CI was reported for each arm. A P value <0.05 was considered significant.</p>		
Clinical significance	- The magnitude of the treatment effect observed was clinically significant.	Yes 1	No 2
Data Contribution Grade (Range 4-8)		4	

Overall S&P Appraisal, Disposition and Weighting


S&P Grade (Range 9-25)	LOE (2) + Suitability (4) + Data Contribution (4) = 10	Disposition and Weighting (select)	Accepted and Pivotal 9-12 Accepted but not Pivotal, 13-21 Excluded, 22-25
---------------------------	---	------------------------------------	--

Relevant S&P Results

Overview	- Among 105 patients who were screened between January 2017 and December 2019, 92 eligible patients (32 women [34.8%]) with a median age of 30 years (IQR: 20-36 years) were randomized equally into the BES and SES groups. The composite of procedural and vascular complications occurred in 10.9% of the BES group and 2.2% of the SES group (odds ratio: 0.18; 95% CI: 0.02-1.62; P = 0.20).
Safety data	<ul style="list-style-type: none"> - The primary composite outcome was observed in five patients (10.9%) in the BES group and one patient (2.2%) in the SES group (OR: 0.18; 95% CI: 0.02-1.62; P = 0.20). Periprocedural complications were observed in three patients (6.5%) in the BES group and one patient (2.2%) in the SES group (OR: 0.31; 95% CI: 0.03-3.18; P = 0.617). - Vascular access complications, consisting of non-flow-limiting femoral artery dissection, femoral artery pseudoaneurysm, and retroperitoneal hemorrhage, occurred with an incidence rate of 1.1% in the overall cohort (two patients [4.3%] in the BES group and no patients in the SES group; P = 0.49). - One patient (1.1%) was complicated by aortic pseudoaneurysm formation, which was subsequently treated with aortic stent-graft implantation (one patient [2.2%] in the BES group, P = 0.31).
Performance data	<ul style="list-style-type: none"> - The procedural success rate was 100%, with no mortality during the 12-month follow-up. The mean catheterization-based postprocedural pressure gradient was 1.4 ± 4.2 mm Hg in the SES group and 1.5 ± 3.2 mm Hg in the BES group (P = 0.52). - Aortic recoarctation was confirmed by cardiac catheterization (pressure gradient >20 mm Hg) in five patients (5.4%) (three patients [6.5%] in the BES group and two patients [4.3%] in the SES group; OR: 0.65; 95% CI: 0.10-4.09; P = 0.64).
Benefits/claims data	- At the one-year follow-up, the median number of antihypertensive medications had dropped from two (IQR: 1-3) to one (IQR: 0-2) in the study population (P < 0.001).
Strengths	- Randomized clinical trial
Weaknesses/ Potential bias	- Limitations: Author-identified limitations include the study may be underpowered, trial was set for 100 patients but reduced to 92 due to shortage of stents, some components of procedural complications were not specifically mentioned in the clinical trial registration website, the severity and clinical impact of vascular access complications varies and may not be comparable, and one year is insufficient for the evaluation of long-term complications.

State of the Art Appraisal

Medical condition	Alternatives	Risk/benefit	Side-effects	Equivalence	Surrogate endpoints
-------------------	--------------	--------------	--------------	-------------	---------------------




NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

	Yes 1	No 2	Yes 1	No 2	Yes 1	No 2	Yes 1	No 2	Yes 1	No 2	Yes 1	No 2										
Overall SOA Appraisal and Disposition																						
SOA Grade (Range 6-12)		10				Disposition (select)				Accepted, < 12 Excluded, 12												
Relevant SOA Results																						
SOA data		- The balloon-expandable stent (BES) was the first and is now the most commonly used aortic stent in CoA treatment with promising results. Nevertheless, the small but considerable risk of stent migration, the unknown true rate of postinterventional aneurysm formation, and the higher needed learning curve for stent application are cited as the most important drawbacks.																				
Comments		- Authors reported that they have no relationships relevant to the contents of this paper to disclose.																				
Safety & Performance Appraisal																						
Level of Evidence		Study Method/Design				Question Applied				Oxford LOE 2011												
		Three-year follow-up to open-label, parallel-group, blinded endpoint randomized pilot clinical trial reported in (18)				To compare the safety and efficacy of the balloon-expandable stent (BES) and the self-expandable stent (SES) in the endovascular treatment of CoA.				<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10%;">1</td> <td style="width: 10%;">2</td> <td style="width: 10%;">3</td> <td style="width: 10%;">4</td> <td style="width: 10%;">5</td> </tr> <tr> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> </table>			1	2	3	4	5					
1	2	3	4	5																		
Suitability		Relevant Data								Grading												
Device		- Uncovered CP BES (n=35 patients) - Uncovered nitinol SES (n=36 patients)								<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10%;">D1</td> <td style="width: 10%;">D2</td> <td style="width: 10%;">D3</td> </tr> <tr> <td></td> <td></td> <td></td> </tr> </table>			D1	D2	D3							
D1	D2	D3																				
Application		- Adult patients with de novo native CoA								<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10%;">A1</td> <td style="width: 10%;">A2</td> <td style="width: 10%;">A3</td> </tr> <tr> <td></td> <td></td> <td></td> </tr> </table>			A1	A2	A3							
A1	A2	A3																				
Patient		- 71 of 92 patients randomized in initial study participated in the three-year structural follow-up (two passed away (one COVID-19 infection, one car accident), three withdrew from study and 16 declined to participate in follow-up) - 25 women (32.2%) with a median age of 30 years (IQR: 20-35 years)								<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10%;">P1</td> <td style="width: 10%;">P2</td> <td style="width: 10%;">P3</td> </tr> <tr> <td></td> <td></td> <td></td> </tr> </table>			P1	P2	P3							
P1	P2	P3																				
Report		- Report suitable for review; the terms thoracic aortic aneurysmal formation and aortic wall injury applied interchangeably.								<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10%;">R1</td> <td style="width: 10%;">R2</td> <td style="width: 10%;">R3</td> </tr> <tr> <td></td> <td></td> <td></td> </tr> </table>			R1	R2	R3							
R1	R2	R3																				
Suitability Grade (Range 4-12)										5												
Data Contribution		Relevant Data								Grading												
Outcomes/Endpoints		- The main outcomes assessed were the three-year rates of recoarctation, aortic injuries, and residual hypertension.								<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10%;">Yes 1</td> <td style="width: 10%;">No 2</td> </tr> <tr> <td></td> <td></td> </tr> </table>			Yes 1	No 2								
Yes 1	No 2																					
Follow-up		- Three-year structural follow-up								<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10%;">Yes 1</td> <td style="width: 10%;">No 2</td> </tr> <tr> <td></td> <td></td> </tr> </table>			Yes 1	No 2								
Yes 1	No 2																					
Statistical analysis		- Data are presented as n (%) or median (IQR). P-values reported for significance.								<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10%;">Yes 1</td> <td style="width: 10%;">No 2</td> </tr> <tr> <td></td> <td></td> </tr> </table>			Yes 1	No 2								
Yes 1	No 2																					
Clinical significance		- The magnitude of the treatment effect observed was clinically significant.								<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10%;">Yes 1</td> <td style="width: 10%;">No 2</td> </tr> <tr> <td></td> <td></td> </tr> </table>			Yes 1	No 2								
Yes 1	No 2																					
Data Contribution Grade (Range 4-8)										4												
Overall S&P Appraisal, Disposition and Weighting																						
S&P Grade		LOE (2) + Suitability (5) +				Disposition and Weighting (select)				Accepted and Pivotal 9-12												



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

(Range 9-25)	Data Contribution (4) = 11		Accepted but not Pivotal, 13-21 Excluded, 22-25
--------------	----------------------------	--	--


Relevant S&P Results

Overview	<ul style="list-style-type: none"> Previously, we reported the one-year results of a randomised controlled trial comparing BES and SES in patients with de novo native CoA. (18) Herein, we have summarised the three-year follow-up results (IRCT20181022041406N3). Of 92 patients initially randomised, 71 patients (25 women [32.2%]), with a median age of 30 years (interquartile range 20-35), participated in the three-year structural follow-up (two patients passed away [one COVID-19 infection and one car accident] and the others did not participate in the follow-up).
Safety data	<ul style="list-style-type: none"> Aortic wall injuries were detected in six patients (8.5%), all treated conservatively with no further endovascular/surgical intervention needed.
Performance data	<ul style="list-style-type: none"> No new recoarctation was detected between the one- and three-year follow-up, and only five patients (with recoarctation previously detected during the first year of follow-up) were identified as having recoarctation. Among those patients, two cases, both initially randomised into the BES group and treated for recoarctation during the first year, needed reballoning due to significant restenosis during the three-year follow-up. We followed up 77.1% (71 of 92) of our randomised population with the structural imaging protocol, and recoarctation occurred in 7.0% of the population with no new cases between the one- and three-year follow-up periods. This finding is in contrast with the major investigations focusing on long-term outcomes, in which a higher rate (~20%) of reintervention has been reported. The inclusion of paediatric patients in the mentioned studies might explain the higher rates of reintervention. Recoarctation rates below 10% were reported when limiting their population to adult patients. A total of 42 out of the 71 patients (59.1%) had residual hypertension, detected more frequently in the BES group, with a trend existing towards a higher median number of antihypertensive drugs during the three-year follow-up. Holzer et al and Eriksson et al reported a downward trend in prolonged hypertension prevalence (42% and 34%, respectively) in patients treated endovascularly. The higher incidence of residual hypertension in the current study might result again from their inclusion of a paediatric population and better blood pressure response in this younger population. <p>Holzer R.J., Gauvreau K., McEnaney K., Watanabe H., Ringel R. Long-Term Outcomes of the Coarctation of the Aorta Stent Trials. Circulation: Cardiovascular Interventions 2021 (582-589) Article Number e010308</p> <p>Eriksson P, Pihkala J, Jensen AS, Dohlen G, Liuba P, Wahlander H, et al. Transcatheter Intervention for Coarctation of the Aorta: A Nordic Population-Based Registry With Long-Term Follow-Up. JACC: Cardiovascular Interventions. 2023;16(4):444-53. doi: 10.1016/j.jcin.2022.11.007.</p>
Benefits/claims data	<ul style="list-style-type: none"> In this three-year follow-up, both BES and SES exhibited low rates of recoarctation, aortic wall injuries and remodeling, but still, more than half of the studied population suffered from residual hypertension.
Strengths	<ul style="list-style-type: none"> Three-year follow-up of randomized clinical trial
Weaknesses/ Potential bias	<ul style="list-style-type: none"> Limitations: Author-identified limitations include small sample size, 23% attrition rate, and lacking ambulatory blood pressure monitoring for residual hypertension.

State of the Art

Appraisal

Medical condition		Alternatives		Risk/benefit		Side-effects		Equivalence		Surrogate endpoints	
Yes 1	No 2	Yes 1	No 2	Yes 1	No 2	Yes 1	No 2	Yes 1	No 2	Yes 1	No 2



NuMED

Summary of Safety and Clinical Performance


SSCP – Stent Placement

Overall SOA Appraisal and Disposition

SOA Grade (Range 6-12)	10	Disposition (select)	Accepted, < 12 Excluded, 12
---------------------------	----	----------------------	--------------------------------

Relevant SOA Results

SOA data	<ul style="list-style-type: none">- Thoracic aortic aneurysmal formation (aortic wall injury) at three years, all treated conservatively with no further endovascular/surgical therapies<ul style="list-style-type: none">o BES: 4/35 (11.4%)o SES: 2/36 (5.6%)
Comments	<ul style="list-style-type: none">- Funding: Study was financially supported by Rajaie Cardiovascular, Medical and Research Center.- Conflict of interest: Authors reported that they have no conflicts of interest to declare.



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

An overall summary of the clinical performance and safety:

A comprehensive, systematic, and critical evaluation of the pertinent clinical data and pre-clinical study data in relation to the BIB Stent Placement catheter has been carried out and documented in the clinical evaluation report. Based on the results of the clinical evaluation report, it is considered that:

- a) Conformity with relevant general safety and performance requirements set out in MDR Annex I under the normal conditions of the intended use of the device has been confirmed.
- b) Undesirable side-effects and acceptability of the benefit-risk ratio have been evaluated and are acceptable according to the current knowledge/the state of the art in the medical fields concerned and according to available medical alternatives.
- c) The information materials supplied by NuMED, and the risk reduction measures are adequate taking into account the intended purpose of the device.
- d) Usability aspects have been adequately considered and the BIB, including the IFUs, is suitable for the intended users.
- e) The claims foreseen in the information materials provided with the CER are adequate taking into account the intended purpose of the device.
- f) The information materials supplied and the RM documentation for the device under evaluation are consistent with the clinical data and pre-clinical study data presented in the CER and with the current knowledge/state of the art.

Overall, it is concluded that the risks associated with the use of the BIB Stent Placement catheter are acceptable when weighed against the benefits to the patient and are compatible with a high level of protection of health and safety, taking into account the generally acknowledged state of the art; that the intended clinical performances are achieved by the device; and that known and foreseeable risks and undesirable side-effects are considered acceptable when weighed against the benefits from performance achieved by the device.

Ongoing planned post-market clinical follow-up:

The BIB Stent Placement Catheter has been commercialized since December 2003 in the EU. Since then, the device is likely to have been used in a variety of patients and populations. A PMCF Study is not warranted at this time due to the fact that long-term safety and clinical performance has been established via device use and ample clinical experience. Continued post-market surveillance activities will provide sufficient data to adequately address clinical risks, and detect emerging risks on the basis of evidence.


The addition of the 26mm, 28mm, and 30mm size in 2014 warranted a PMCF Study just for those sizes. The study investigated procedural and device complications and compared that with currently available data from the use of the smaller BIB Stent Placement Catheter sizes. Based on the findings from the PMCF study, NuMED determined the BIB Stent Placement Catheter (larger sizes) to be safe and effective when used for the approved indication. No changes were required to the risk analysis as there were no new risks identified, and no changes were required to the instructions for use based on the results of the PMCF Study.

6. Possible diagnostic or therapeutic alternatives

Alternative therapies to balloon dilatation / stenting include balloon angioplasty without a stent and surgical intervention.

7. Suggested profile and training for users

The device is intended to be used by a cardiac surgeon and/or interventionalist.



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

8. Reference to any harmonised standards and CS applied


There are no Common Specifications for this type of device.

The following harmonised standards are followed for this device:

- EN ISO 10993-10: 2023 – Biological Evaluation of Medical Devices – Part 10: Tests for Skin Sensitization
- EN ISO 10993-18: 2020 – Biological Evaluation of Medical Devices – Part 18: Chemical characterization of medical device materials within a risk management process
- EN ISO 10993-23: 2021 – Biological Evaluation of Medical Devices – Part 23: Tests for Irritation
- EN ISO 11135: 2014 / A1:2019 – Sterilization of health-care products – Ethylene oxide – Requirements for the development, validation and routine control of a sterilization process for medical devices.
- BS EN ISO 11607-1: 2020 +A1: 2023 – Packaging for Terminally Sterilized Medical Devices – Part 1: Requirements for materials, sterile barriers systems and packaging systems
- BS EN ISO 11607-2: 2020 +A1: 2023 – Packaging for Terminally Sterilized Medical Devices – Part 2: Validation requirements for forming, sealing and assembly processes
- EN ISO 11737-1: 2018 / A1:2021 – Sterilization of medical devices – Microbiological methods – Part 1: Determination of a population of microorganisms on products
- EN ISO 13485: 2016 / A11:2021 – Medical devices – Quality management systems – Requirements for regulatory purposes
- EN ISO 14971: 2019 / A11:2021 – Medical Devices – Application of Risk Management to Medical Devices
- EN ISO 15223-1: 2021 – Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements

9. References

1. Sohrabi B, Jamshidi P, Yaghoubi A, Habibzadeh A, Hashemi-Aghdam Y, Moin A, Kazemi B, Ghaffari S, Abdolazadeh BM, Mahmoodi K, Comparison between covered and bare Cheatham-Platinum stents for endovascular treatment of patients with native post-ductal aortic coarctation: immediate and intermediate-term results. *JACC. Cardiovascular interventions* **7**(4), 416-423 (2014).
2. Vanagt WY, Cools B, Boshoff DE, Frerich S, Heying R, Troost E, Louw J, Eyskens B, Budts W, Gewillig M, Use of covered Cheatham-Platinum stents in congenital heart disease. *International Journal of Cardiology* **175**, 102-107 (2014).
3. Alcibar J, Blanco R, Fernandez L, Arriola J, Garcia K, Pena N, Inguanzo R, Voces R, Castellanos E, Montes PM, Elective implantation of covered stents for coarctation and recoarctation in adolescents and adults. *Revista espanola de cardiologia (English ed.)* **66**, 443-449 (2013).
4. Chang ZP, Jiang SL, Xu ZY, Zhang GJ, Huang LJ, Zhao SH, Ling J, Zheng H, Jin JL, Wu WH, Hu HB, Li SG, Yu JH, Yan CW, Use of covered Cheatham-Platinum stent as the primary modality in the treatment for native coarctation of the aorta. *Chinese medical journal* **125**, 1005-1009 (2012).
5. Erdem A, Akdeniz C, Sarıtaş T, Erol N, Demir F, Karaci AR, Yalçın Y, Celebi A, Cheatham-Platinum stent for native and recurrent aortic coarctation in children and adults: immediate and early follow-up results. *Anadolu Kardiyol Derg.* **Aug;11(5)**, 441-449 (2011).
6. Butera G, Heles M, MacDonald ST, Carminati M, Aortic coarctation complicated by wall aneurysm: the role of covered stents. *Catheterization and cardiovascular interventions : official journal of the Society for Cardiac Angiography & Interventions* **78**, 926-932 (2011).
7. Moltzer E, Roos-Hesselink JW, Yap SC, Cuypers JA, Bogers AJ, de Jaegere PP, Witsenburg M, Endovascular stenting for aortic (re)coarctation in adults. *Netherlands heart journal : monthly journal of the Netherlands Society of Cardiology and the Netherlands Heart Foundation* **18**, 430-436 (2010).
8. Kische S, Schneider H, Akin I, Ortak J, Rehders TC, Chatterjee T, Nienaber CA, Ince H. Technique of interventional repair in adult aortic coarctation. *J Vasc Surg* **Jun;51(6)**, 1550-1559 (2010).
9. Agnoletti G, Marini D, Ou P, Vandrell MC, Boudjemline Y, Bonnet D, Cheatham platinum (CP) and Palmaz stents for cardiac and vascular lesions treatment in patients with congenital heart disease. *EuroIntervention* **4**, 620-625 (2009).
10. Peters B, Ewert P, Berger F, The role of stents in the treatment of congenital heart disease: Current status and future perspectives. *Ann Pediatr Cardiol* **Jan;2(1)**, 3-23 (2009).
11. Tzifa A, Ewert P, Brzezinska-Rajszys G, Peters B, Zubrzycka M, Rosenthal E, Berger F, Qureshi SA, Covered Cheatham-platinum stents for aortic coarctation: early and intermediate-term results. *J Am Coll Cardiol* **47**, 1457-1463 (2006).
12. Cheatham JP, Stenting of coarctation of the aorta. *Catheterization and cardiovascular interventions : official journal of the Society for Cardiac Angiography & Interventions* **54**, 112-125 (2001a).
13. Cheatham JP, Improved stents for pediatric applications. *Progress in Pediatric Cardiology*, **14**, 95-115 (2001b).



NuMED

Summary of Safety and Clinical Performance

SSCP – Stent Placement

14. Meadows J, Minahan M, McElhinney DB, McEnaney K, Ringel R, Intermediate Outcomes in the Prospective, Multicenter Coarctation of the Aorta Stent Trial (COAST). *Circulation* **131**, 1656-1664 (2015).
15. Bairam AR, Ali AH, Al Ghizzi HJ. Immediate and intermediate term results of Cheatham platinum Stenting for native coarctation of aorta in adults. *Revista Latinoamericana de Hipertension*. 2021;15(3):192-9. doi: 10.5281/zenodo.5651255.
16. Schleiger A, Al Darwish N, Meyer M, Kramer P, Berger F, Nordmeyer J. Long-term follow-up after endovascular treatment of aortic coarctation with bare and covered Cheatham platinum stents. *Catheterization and Cardiovascular Interventions*. 2023;102(4):672-82. doi: 10.1002/ccd.30793.
17. Vargas-Acevedo C, Mejia E, Kagiya Y, Soszyn N, Zablah JE, Morgan GJ. Use of a mechanical crimper for mounting covered stents in aortic coarctation. *Progress in Pediatric Cardiology*. 2024;74. doi: 10.1016/j.ppedcard.2024.101730.
18. Sadeghipour P, Mohebbi B, Firouzi A, Khajali Z, Saedi S, Shafe O, et al. Balloon-Expandable Cheatham-Platinum Stents Versus Self-Expandable Nitinol Stents in Coarctation of Aorta: A Randomized Controlled Trial. *JACC: Cardiovascular Interventions*. 2022;15(3):308-17. doi: 10.1016/j.jcin.2021.11.025.
19. Sadeghipour P, Pouraliakbar HR, Farrashi M, Habibi Khorasani S, Mohebbi B, Iranian M, et al. Balloon-expandable versus self-expanding stents in native coarctation of the aorta: three-year results of a randomised controlled trial. *EuroIntervention*. 2024;20(9):613-5. doi: 10.4244/eij-d-23-00846. PubMed PMID: 38726718; PubMed Central PMCID: PMC11067719.

10. Revision History

SSCP revision number	Date Issued	Change Description	Revision validated by Notified Body
00	21 June 2022	Initial implementation	<input type="checkbox"/> Yes Validation Language: English <input checked="" type="checkbox"/> No
01	07 July 2023	Updated sections 5, 8, and 9 for CER Update.	<input type="checkbox"/> Yes Validation Language: English <input checked="" type="checkbox"/> No
02	01 August 2025	Updated section 1 (EU Rep info), Section 2 limitations, Section 3 previous generations, variants and accessories, Section 5 for updated clinical literature, Section 7 users, Section 8 for harmonized standards and Section 9 for references.	<input type="checkbox"/> Yes Validation Language: English <input checked="" type="checkbox"/> No