MRI Safety Information

Nonclinical testing and modeling has demonstrated that the Stent is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3 T
- Maximum spatial gradient magnetic field of 2500 gauss/cm (25 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg for 15-minutes of scanning (Normal Operating Mode)

Based on nonclinical testing and modeling, under the scan conditions defined above, the Stent is expected to produce a maximum in vivo temperature rise of less than 2°C after 15 minutes of continuous scanning.

MR image quality may be compromised if the area of interest is in the same area, or relatively close to the position of the device. In nonclinical testing, the image artifact caused by the device extends approximately 3 mm from the Stent when imaged with a spin echo pulse sequence and 6 mm when imaged with a gradient echo pulse sequence and a 3 T MRI System. The lumen of the device was obscured.

The presence of other implants or medical circumstances of the patient may require lower limits on some or all of the above parameters.

The covered stent is composed of heat treated 90% platinum / 10% iridium wire that is arranged in a "zig" pattern, laser welded at each joint, and over brazed with 24K gold. The covered stent has an ePTFE covering attached to the stent framework with a cyanoacrylate adhesive.

Device Lifetime: minimum 10 years.

Covered Stent Implant Card

NuMED, Inc.
2880 Main Street
Hopkinton, NY 12965
1-877-NuMEDCP(686-3327) or
(315) 328-4491

AFFIX LABEL WITH UDI INFORMATION

INSTRUCTIONS

Physician: Please fill out (1) card and return to NuMED, Inc., 2880 Main St., Hopkinton, NY 12965, USA; ATTN: RA Mgr. Please fill out (1) card and give to patient.

Patient: Please carry this card at all times and show it to any medical personnel who may be treating you. The purpose of this card is to enable NuMED, Inc. to notify the patient of new information regarding the safety and effectiveness or performance of the implant, and any required corrective action. Please advise NuMED of any change of address.

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