
CE
1639

MULTI-TRACK™

ANGIOGRAPHIC CATHETER

Instructions for Use.....	2
Mode d'emploi	3
Istruzioni per l'uso.....	5
Gebrauchsanweisung	6
Instrucciones de utilización	8
Bruksanvisning	9
Brugsanvisning	11
Gebruiksaanwijzing.....	12
Instruções de utilização	14
Kullanım Talimatları	15
Οδηγίες χρήσης.....	17

CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician.

Instructions for Use

INDICATIONS

Recommended for use in catheterization for angiography of cardiovascular vessels and / or chambers. It can be used for injection of contrast medium and pressure measurement in any chamber or vessel.

CATHETER CHARACTERISTICS

The MULTI-TRACK™ angiographic catheter is a single lumen catheter with side holes at the tip and a 1cm section at the distal tip for catheter tracking. The tubing of the catheter is radiopaque for proper visualization under fluoroscopy.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas. Sterile and non-pyrogenic if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

WARNING

- This device is intended for single use only. Do not resterilize and/or reuse it, as this can potentially result in compromised device performance and increased risk of cross-contamination.
- Do not advance the MULTI-TRACK™ angiographic independent from the guidewire.
- Do not advance the guidewire or any other component if resistance is met, without first determining the cause and taking remedial action.
- The Rated Burst Pressure for this device is 1000 psi. Do Not Exceed 1000 psi.

PRECAUTIONS

- The MULTI-TRACK™ angiographic catheter procedures should be conducted under fluoroscopic/MRI guidance with appropriate equipment.
- Guidewires are delicate instruments. Care should be exercised while handling to help prevent the possibility of breakage.
- Careful attention must be paid to the maintenance of tight catheter connections and aspiration before proceeding to avoid air introduction into the system.
- Under no circumstances should any portion of the catheter system be advanced against resistance. The cause of the resistance should be identified with fluoroscopy/MRI and action taken to remedy the problem.
- Proper functioning of the catheter depends on its integrity. Care should be used when handling the catheter. Damage may result from kinking, stretching, or forceful wiping of the catheter.
- Do not use excessive force on the catheter. During an angiographic procedure, excessive force is determined to be the force required to bend a standard 0.035" guidewire.
- The flow rates have been calculated as follows at 1000 psi:

Multi-Track™ Angiographic Catheter Flow Chart

Catheter Size (Fr)	Catheter Length (cm)	Flow Rate (ml/sec)
3	80	5.5
3	100	4
4	80	13
4	100	11
5	80	22
5	100	20
6	100	25

MRI Safety Information

Non-clinical testing has demonstrated that NuMED's Catheter Systems are MR Conditional. The following labeling is specific to NuMED's catheters and does not apply to guidewires or other accessory devices used in combination with NuMED's Catheter Systems. Users should consult the instructions for use and MRI safety information of other accessory devices used with NuMED's catheters. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5-Tesla (1.5 T) or 3-Tesla (3 T).
- Maximum spatial field gradient of 1,900 G/cm (19 T/m).
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4.0 W/kg (First Level Operating Mode).

RF Heating

Under the scan conditions defined above, NuMED's Catheter Systems are expected to produce a maximum temperature rise of less than 5.7°C after 15 minutes of continuous scanning. Caution: The RF heating behavior does not scale with static field strength. Devices that do not exhibit detectable heating at one field strength may exhibit high values of localized heating at another field strength.

MR Artifact

In non-clinical testing, the image artifact caused by the catheter extend radially approximately 3 mm from the catheter when imaged in a 3 T MRI system.

POTENTIAL COMPLICATIONS

Potential complications & adverse effects associated with device use and indication include:

- (1) infection,
- (2) air embolism, and
- (3) hematoma formation.

Any serious incident that has occurred in relation to the device should be reported to NuMED & the Competent Authority in the country of use.

PREPARATORY MEASURES

1. Remove the catheter from its packaging and carefully examine for bends or kinks. Carefully flush the MULTI-TRACK™ angiographic catheter with saline.

INSTRUCTIONS

1. Enter the vessel percutaneously using the standard Seldinger technique over the appropriate guidewire for the size catheter being used.
2. Position an appropriate guidewire for the catheter being used into the desired site. Connect the distal portion of the MULTI-TRACK™ angiographic catheter to the guidewire and advance the catheter to the entrance into the body. Flush the catheter in order to remove all air bubbles from its lumen.
3. Hold the guidewire and gently advance the MULTI-TRACK™ angiographic catheter into the skin. If resistance is encountered, put an appropriate guidewire into the lumen of the MULTI-TRACK™ angiographic catheter in order to stiffen the shaft while pushing.
4. Once the catheter is in the vessel, flush again according to standard practice.
5. Advance the catheter into the desired position either by simply pushing on its shaft or with the aid of a guidewire in its lumen.
6. Perform pressure measurements and/or angiographies according to the standard practice or hospital protocols.
7. After terminating the use of the MULTI-TRACK™ angiographic catheters, they are sequentially removed and another procedure can be performed on the guidewire which was positioned in the beginning.
8. After removal of all devices from the vascular access apply pressure to the insertion site according to standard practice or hospital protocol for percutaneous vascular procedures.
9. Dispose of device after use according to standard hospital protocol for biohazardous devices.

WARNING: NuMED catheters are placed in the extremely hostile environment of the human body. Catheters may fail to function for a variety of causes including, but not limited to, medical complications or failure of catheters by breakage. In addition, despite the exercise of all due care in design, component selection, manufacture and testing prior to sale, catheters may be easily damaged before, during, or after insertion by improper handling or other intervening acts. Consequently, no representation or warranty is made that failure or cessation of function of catheters will not occur or that the body will not react adversely to the placement of catheters or that medical complications will not follow the use of catheters.

NuMED cannot warrant or guarantee NuMED accessories because the structure of the accessories may be damaged by improper handling before or during use. Therefore, no representations or warranties are made concerning them.

Warranty and Limitations

Catheters and accessories are sold in an 'as is' condition. The entire risk as to the quality and performance of the catheter is with the buyer. NuMED disclaims all warranties, expressed or implied, with respect to catheters and accessories, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose. NuMED shall not be liable to any person for any medical expenses or any direct or consequential damages resulting from the use of any catheter or accessory or caused by any defect, failure, or malfunction of any catheter or accessory, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind NuMED to any representation or warranty with respect to catheters and accessories.

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
REFERENCES

- Youkey Jr, Clagett GP, Rich NM, Jaffin JH, Cohen AJ, Brigham RB, Orecchia PM, Saalander JM, Vascular Trauma secondary to diagnostic and therapeutic procedures 1974 through 1982 A comparative review. *AM J Surg* 146(6): 788-791, 1983
- Mahnke PF, Pathology of vein, artery and heart catheterization Studies of autopsies. *Z Gesamte Inn Med.* 41 (13): 378-382, 1986.
- Ducatman BS, McMichan JC, Edwards WD Catheter induced lesions of the right side of the heart. A one-year prospective study of 141 autopsies. *JAMA* 253 (6):791-795, 1985.
- Lange HW, Galliani CA, Edwards JE. Local complications associated with indwelling Swan-Ganz catheters: autopsy study of 36 cases. *Am J Cardiol.* 32 (18): 1108-1, 1983.
- McCullum Ch, Mavor E, Brachial artery injury after cardiac catheterization. *J Vasc. Surg.* 4 (4) 355-359, 1986. Champault G, Totally implantable catheters for cancer chemotherapy: French experience on 325 cases. *Cancer Drug Deliv* 3 (2): 131, 1986.
- Wolte BM, Ryder MA, Nishikawa RA, Halsted CH, Schmidt BF, Complications of parenteral nutrition. *Am J Surg* 152 (1) 93-9, 1986.
- Jensen Sr, Crummy AB, McDermott JC, Starck E, Voegeli DR, DSA evaluation of indwelling central venous catheter. *Cardiovasc. Interent Radiol* 8 (4): 216-2177. M1985.
- Kido DK, Paulin S, Alenghat JA, Waternaux C, Riley WD, Thrombogenicity of heparin and nonheparin coated catheters: clinical trial. *AJR* 139 (5): 957-961, 1982.
- Coronary spasm producing coronary thrombosis. *N Engl J Med* 309 (26): 1648, 1983.
- Bonhoeffer, Philipp; Piechaud, Jean-Francois; The Multi-Track angiography catheter: a new tool for complex catheterisation in congenital heart disease. (76); 173-177, (1996).

	Do Not Reuse Ne pas réutiliser Non riutilizzare Nicht wiederverwenden No reutilizar	Φάρ ει återanvändas Må ikke genbruges Niet hergebruiken NÃO Reutilizar Tekrar Kullanmayın Μην επαναχρησιμοποιείτε		Catalogue Number Número de catálogo Katalognummer Numero di catalogo Καταλόγος αριθμός Número de catálogo	Katalognummer Καταλόγος αριθμός Καταλόγος αριθμός Número de Catálogo Καταλόγος αριθμός Αριθμός καταλόγου
	Batch Code Code du lot Codice del lotto Chargenode Código del lote	Σάνδινγκκωδ Partikode Partijcode Código do Lote Parti Kodu Κωδικός παρτίδας		Use By À utiliser avant Utilizzare entro Verwendbar bis Usar antes de	Αναvänd senast Anvendes inden Gebruiken voor Usar Até Son Kullanma Tarihi Ημερομηνία λήξης
	Do not resterilize Ne pas résteriliser Non risterilizzare Nicht resterilisieren No reesterilizar			Får inte omsteriliseras Må ikke resteriliseres Niet hersteriliseren NÃO reesterilizar Yeniden sterilize etmeyin Μην επαναστεριώνεται	
	Attention: See Instructions For Use Attention : Consulter le mode d'emploi Attenzione – Consultare le istruzioni per l'uso Achtung: Siehe Gebrauchsanweisung Atención: consulte instrucciones de uso			OBS! Se bruksanvisningen Bemærk: Se brugsanvisningen Let op: zie gebruiksaanwijzing Atenção: Consulte as instruções de utilização Dikkat: Kullanım Talimatlarına Bakın Προσοχή: Δείτε τις οδηγίες χρήσης	
	Keep Away From Sunlight Garder à l'abri de la lumière du soleil Proteggere dai raggi solari Sonneneinstrahlung vermeiden Mantener alejado de la luz solar			Exponera inte för direkt solljus Holdes væk fra sollys Uit direct zonlicht houden Manter Alastado da Luz Solar Güneş Işığından Uzak Tutun Διατηρείτε μακριά από ηλιακό φως	
	Temperature Limitation Limite de température Limitazione di temperatura Temperaturbegrenzung Limitación de la temperatura			Temperaturbegränsning Temperaturbegrænsning Temperaturlæperking Limite de Temperatura Isi Sınırlaması Περιορισμός θερμοκρασίας	
	Manufacturer / Date of Manufacture Fabricant / Date de fabrication Produttore / Data di fabbricazione Hersteller / Herstellungsdatum Fabricante / Fecha de fabricación			Tilværkare / Tilværningsdatum Producent / fremstillingsdato Fabrikant / Productiedatum Fabricante / Data de Fabrico Üretici / Üretim Tarihi Κατασκευαστής / Ημερομηνία κατασκευής	
	Sterilized Using Ethylene Oxide Stérilisé à l'oxyde d'éthylène Sterilizzato con ossido di etilene Sterilisiert mit Ethylenoxid Esterilizado con óxido de etilén			Steriliserad med etylenoxid Steriliseret med ethylenoxid Gesteriliseerd met ethyleenoxide Esterilizado Usando Óxido de Etileno Etilen Oksit Kullanılarak Sterilize Edilmiştir Αποστειρώθηκε με χρήση οξείδιου του αιθυλενίου	
	Authorised Representative in the European Community Représentant autorisé dans la Communauté européenne Rappresentante autorizzato nella Comunità Europea Autorisierter Vertreter in der Europäischen Gemeinschaft Representante autorizado en la Comunidad Europea			Auktoriserad representant inom EU Autoriseret repræsentant i Det Europæiske Fællesskab Gevoimachtigd vertegenwoordiger in de Europese Gemeenschap Representante Autorizado na Comunidade Europeia Avruppa Topluluğunda Yetkili Temsilci Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα	
	Do not use if package is damaged Ne pas utiliser si l'emballage est endommagé Non utilizzare se l'imballaggio non è integro Nicht verwenden, wenn die Verpackung beschädigt ist No utilizar si el envase está dañado			Får inte användas om förpackningen är skadad Må ikke anvendes, hvis emballagen er beskadiget Niet gebruiken als de verpakking beschadigd is NÃO utilizar se a embalagem estiver danificada Ambalaj hasarlıysa kullanmayın Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά	
	MR Conditional Utilisable(s) en environnement IRM sous certaines conditions A compatibilità RM condizionata Bedingt MR-sicher Producto seguro bajo ciertas condiciones de la resonancia magnética			Villkorligt godkänd för MR Betinget MR-sikker MRI-veilig onder bepaalde voorwaarden MR Conditional (pode ser submetido a exames de RM sob determinadas condições) MR Koşullu Ασφαλές για μαγνητική τομογραφία υπό προϋποθέσεις	
	Medical Device Dispositif Médical Dispositivo medico Medizinprodukt Producto sanitario			Medicinteknisk produkt Medicinsk udstyr Medisch hulpmiddel Dispositivo médico Tıbbi Cihaz Ιατροτεχνολογικό προϊόν	

IFU-310CE-16

NuMED, Inc.
2880 Main Street
Hopkinton, NY
USA 12965
TEL: (315) 328-4491
FAX: (315) 328-4941
email: info@numedusa.com
www.numedforchildren.com

European Representative:
G. van Wageningen B.V.
Hallenweg 40, 5683 CT Best
The Netherlands
TEL: +31-499-377388
FAX: +31-499-377456
email: info@heartmedical.nl

RM0344-171

20 November 2020