
CE  
1639

# MITRAL DILATATION KIT™

## BONHOEFFER MULTI-TRACK SYSTEM

Instructions for Use .....	2
Mode d'emploi.....	5
Istruzioni per l'uso .....	8
Gebrauchsanweisung.....	11
Instrucciones de utilización.....	14
Bruksanvisning.....	17
Brugsanvisning .....	20
Gebbruksaanwijzing .....	23
Instruções de utilização.....	26
Οδηγίες χρήσης .....	29

**CAUTION:** Federal (USA) Law restricts this device to sale by or on the order of a physician.

## Instructions for Use

**INDICATIONS / INCLUSION**

The Bonhoeffer MULTI-TRACK™ Mitral Dilatation Kit is indicated for Percutaneous Transvenous Mitral Commissurotomy in patients with hemodynamically significant mitral valvular stenosis resulting primarily from commissural fusion of the mitral valve cusps.

**DESCRIPTION**

The Bonhoeffer MULTI-TRACK™ Mitral Kit features two balloon dilatation catheters, one 5F MULTI-TRACK™ angiographic catheter, one 14F dilator, and one preformed .035" guidewire. The MULTI-TRACK™ balloon catheter is threaded over a guidewire with the most distal 1cm beyond the balloon. This will allow a rapid exchange balloon catheter to be threaded over the guidewire in the usual way to allow the balloons to be side by side in the mitral valve.

The shafts of the balloon catheters are marked with the balloon size and serve as the balloon inflation and deflation lumen. The balloons are marked with two radiopaque markers beneath the shoulders. Each balloon inflates to the stated diameter and length at a specific pressure. The balloon size is  $\pm 10\%$  at the Rated Burst Pressure (RBP). The RBP is different for each size. Please check the package label for the RBP. It is important that the balloon not be inflated beyond the RBP. Each catheter is intended for a single application and should not be reused.

**HOW SUPPLIED**

Supplied sterilized by ethylene oxide gas. Sterile and non-pyrogenic if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

**CONTRAINDICATIONS / EXCLUSION**

The patient's medical conditions can affect successful use of these catheters. Standard risks of trans-septal catheterization and the risk of erroneous positioning of the guidewire as well as erroneous use of the balloon catheter have to be considered. Further, the dilatation of a heavily calcified valve can favor balloon rupture even under the RBP.

- Patients with mild stenosis and no functional impairment.
- Patients who require open heart surgery for the treatment of other valves.
- Patients with significant (>++) mitral regurgitation.
- Patients with left atrial thrombus.
- Patients who have an excessively high echo-score [1]

**WARNINGS**

- The use of this kit must be carefully considered in patients with heavy calcifications.
- The inflation diameter of both balloons must be carefully considered in selecting the sizes for a particular patient. The sum of the balloon sizes should not be significantly bigger than the annulus of the mitral valve. [2]
- Caution: Do not exceed the RBP. An inflation device with a pressure gauge is recommended to monitor pressure.
- Use only appropriate balloon inflation medium. Do not use air or gaseous medium to inflate the balloon.
- Do not remove the guidewire from the catheter at any time during the procedure.
- When using a guidewire, potential exists for thrombus formation or emboli, arterial or venous wall damage and/or plaque dislodgment. The physician should be familiar with the literature concerning the complications of angiography.
- Exercise care in handling of the guidewire during a procedure to reduce the possibility of accidental breakage, bending, kinking, or coil separation. Do not use a guidewire that has been damaged. Do not attempt to straighten a wire that has been kinked or bent. Do not advance a kinked guidewire into a balloon catheter to reduce the potential of wire breakage. Follow the recommendation about guidewire placement described by Dr. Bonhoeffer [3].
- This device is intended for single use only. Do not resterilize and/or reuse it, as this can potentially result in compromised device performance and increased risk of cross-contamination.
- These catheters should be used prior to the 'Use Before' date noted on the package label.
- DO NOT EXCEED the RBP stated on the package label.
- This kit is intended for mitral valvuloplasty procedures, and is not intended for coronary angioplasty.

**PRECAUTIONS**

- Proper catheter function depends upon its integrity. Care must be taken to prevent damage during handling and operation from kinking, stretching, wiping with excessive force, guidewire insertion and advancement, or other procedures involving its use.
- Careful attention must be paid to establish and maintain tight connections in order to prevent air introduction into the catheters or fluid leakage from the catheters during its operation.
- Guidewires are delicate instruments. Care should be exercised while handling to help prevent the possibility of breakage.
- To prevent possible tissue damage, care should be taken when manipulating a catheter over a guidewire during the balloon catheter's placement and withdrawal. If resistance is felt during catheter placement, discontinue the procedure and determine the cause of resistance before proceeding. If the cause of resistance cannot be determined, remove the guidewire and catheter as a unit to prevent possible damage and/or complications.
- Dilatation procedure should be conducted under fluoroscopic guidance with appropriate x-ray equipment.
- Dr. Bonhoeffer recommends performing a percutaneous puncture without a sheath. However, when using a sheath if resistance is felt upon removal, then the balloons, guidewire and the sheath should be removed together as a unit, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the balloon catheters and sheath as a unit and withdrawing both together, using a gentle twisting motion combined with traction.
- Before removing the balloon catheters it is very important that the balloons are completely deflated.

**POTENTIAL COMPLICATIONS/ADVERSE EFFECTS**

Potential balloon separation following balloon rupture or abuse and the subsequent need to use a snare or other medical interventional techniques to retrieve the pieces.

**NOTE:** There have been infrequent reports of larger diameter balloons bursting circumferentially, possibly due to a combination of tight focal strictures in large vessels. In any instance of a balloon rupture while in use, it is recommended that a sheath be placed over the ruptured balloon prior to withdrawal through the entry site. This can be accomplished by cutting off the proximal end of the catheter and slipping an appropriately sized sheath over the catheter into the entry site. For specific technique, refer to: [4].

Potential complications & adverse effects associated with device use and indication include:

- Perforation
- Hematoma
- Arrhythmia Development
- Valvular Tearing or Trauma
- Restenosis Development
- Inflammation
- Infection
- Septal Puncture

- Embolic Events
- Cardiac Tamponade
- Hemopericardium
- Mitral Regurgitation
- Tachycardia
- Death

Guidewire complications: Vascular Damage, Thrombus Formation, & guidewire perforation with cardiac tamponade.

Any serious incident that has occurred in relation to the device should be reported to NuMED & the Competent Authority in the country of use.

**INSPECTION AND PREPARATION**

Prior to valvuloplasty, carefully examine all equipment to be used during the procedure, including the catheters, to verify proper function and that the catheter sizes are suitable for the specific procedure for which it is intended. Also, inflate the dilatation catheters to the appropriate RBP and deflate to verify proper function.

1. Remove the balloon protectors and inspect the catheters for bends or kinks prior to use.
2. Prepare a 50% mixture by volume of contrast medium and normal saline.
3. Purge the catheters through lumen thoroughly, observing for leaks.
4. To check inflation/deflation times, use a stopwatch. Start the inflation cycle. Begin timing at the beginning of the inflation cycle. The balloons should appear fully inflated within 5 seconds. Deflate the balloons. The balloons are considered totally deflated when the balloon has collapsed against the catheter shaft and contrast medium is no longer visible. A deflation time of 10 seconds or less is acceptable. Repeat the procedure several times to verify the inflation/deflation time.
5. Point the inflator nozzle downward and aspirate until all air is removed from the balloons, and bubbles no longer appear in the contrast solution.
6. Turn the stopcock off to maintain the vacuum in the balloons.

**INSERTION: Vascular**

1. Enter the vessel percutaneously using the standard Seldinger technique over the appropriate guidewire for the size catheter being used.
2. Advance the catheter across the lesion with fluoroscopic guidance using accepted percutaneous transluminal angioplasty technique. In most patients, the balloon should meet with minimal resistance to insertion. Do not advance the catheter unless the guidewire is in place.
3. Referring to the chart below, inflate the balloons with contrast medium until the desired diameter is achieved or the RBP is reached whichever comes first. Do Not Exceed the RBP.

**DEFLATION AND WITHDRAWAL**

1. Deflate the balloons by drawing a vacuum with an inflator. Note: The greater the vacuum applied and held during withdrawal, the lower the deflated balloon profile.
2. Gently withdraw the catheters. As the balloon exits the vessel, use a smooth, gentle, steady, TRACTION motion. If resistance is felt upon removal, then the balloon and the sheath should be removed together as a unit under fluoroscopic guidance, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the balloon catheters and sheath as a unit and withdrawing both together, using a gentle twisting motion combined with traction.
3. Apply pressure to the insertion site according to standard practice or hospital protocol for percutaneous vascular procedures.
4. Dispose of device after use according to standard hospital protocol for biohazardous devices.

**Multi-Track Balloon & Rapid Exchange Catheter Sizing Chart**

Applied Press.	14.0 (mm)	16.0 (mm)	18.0 (mm)	20.0 (mm)
1.0 ATM	13.40	15.20	16.80	19.40
2.0 ATM	13.70	15.60	17.40	19.90
3.0 ATM	13.90	16.00	18.00	20.60
4.0 ATM	14.20	16.40	<b>18.60</b>	<b>21.10</b>
5.0 ATM	14.40	<b>16.80</b>		
6.0 ATM	<b>14.60</b>			

The figures in bold face represent the balloon diameter @ Rated Burst Pressure.

**FOR ALL NuMED CATHETERS AN INFLATION DEVICE WITH PRESSURE GAUGE SHOULD BE USED.**

**WARNING:** NuMED catheters are placed in the extremely hostile environment of the human body. Catheters may fail to function for a variety of causes including, but not limited to, medical complications or failure of catheters by breakage. In addition, despite the exercise of all due care in design, component selection, manufacture and testing prior to sale, catheters may be easily damaged before, during, or after insertion by improper handling or other intervening acts. Consequently, no representation or warranty is made that failure or cessation of function of catheters will not occur or that the body will not react adversely to the placement of catheters or that medical complications will not follow the use of catheters.

NuMED cannot warrant or guarantee NuMED accessories because the structure of the accessories may be damaged by improper handling before or during use. Therefore, no representations or warranties are made concerning them.

**Warranty and Limitations**

Catheters and accessories are sold in an 'as is' condition. The entire risk as to the quality and performance of the catheter is with the buyer. NuMED disclaims all warranties, expressed or implied, with respect to catheters and accessories, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose. NuMED shall not be liable to any person for any medical expenses or any direct or consequential damages resulting from the use of any catheter or accessory or caused by any defect, failure, or malfunction of any catheter or accessory, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind NuMED to any representation or warranty with respect to catheters and accessories.

## Instructions for Use of the Bonhoeffer MULTI-TRACK™ Mitral Dilatation Kit

Description of the procedure for mitral dilatation:

1. After local anesthesia the right femoral vein is punctured and a standard trans-septal procedure is performed with a 6F Mullins introducer. **(Figure 1)**
2. After trans-septal puncture the patient is heparinized.
3. The dilator of the Mullins sheath is removed and a balloon wedge catheter with end hole, (supplied by doctor) is introduced through the sheath and into the mitral orifice. **(Figure 2)**
4. Advance the guidewire through the balloon wedge catheter into the left ventricle. Full deployment of the floppy part of the guidewire towards the left ventricular outflow is mandatory before removing the Mullins sheath and the balloon catheter. Adequate guide position must be maintained throughout the procedure [3].
5. Hemodynamic appreciation of the trans-mitral pressure gradient is then obtained by simultaneous measurement of the left ventricular and left atrial pressure. This can be obtained through a single venous approach using two MULTI-TRACK™ Angiographic catheters which are advanced over the previously positioned guidewire, the first in the left ventricle and the second in the left atrium. In most instances the catheters readily advance through the atrial septum, if not the catheters can be stiffened by introducing a guidewire into the lumen in order to enhance pushability. If a pigtail catheter has been positioned in the aortic root for trans-septal puncture this catheter can be advanced into the left ventricle and left atrial pressure can be measured either through the Mullins sheath or by using one single MULTI-TRACK™ Angio catheter advanced into the left atrium.
6. After the assessment of the trans-mitral pressure gradient the MULTI-TRACK™ Angiographic catheters are removed and the skin and atrial septum are dilated with a 14F dilator. **(Figure 3)**
7. The choice of the balloons is made in a manner that the sum of their diameters is approximately 90-100% of the mitral annulus measured by echocardiography [3].
8. The MULTI-TRACK™ balloon catheter tip, which will run on the guidewire is then connected with the guidewire. The catheter is then advanced as described above for the MULTI-TRACK™ angiographic catheter. The catheter is pushed through the mitral orifice. The radiopaque bands mark the shoulders of the dilatation balloons. Avoid kinking of the dilatation balloons. Avoid kinking of the catheter shaft at cutaneous level or during crossing of the atrial septum.
9. When the first balloon is in place a dilatation is performed in order to give some relief to patients with severe mitral stenosis preventing further reduction of the mitral valve area by the presence of this catheter. **(Figure 4)**
10. The second balloon catheter, which is a rapid exchange dilatation catheter, is advanced over the guidewire and lined up with the first balloon.
11. A manual simultaneous inflation of both balloons is then performed at a pressure sufficient to make the "waist" disappear. **(Figure 5)**
12. After the dilatation the balloons are deflated and removed sequentially leaving the guidewire in place.
13. Simultaneous pressure measurements with the help of one or two MULTI-TRACK™ Angiographic catheters are then repeated as described above, in order to evaluate the result.
14. If the result is judged not to be sufficient it is possible to repeat the procedure with either a larger MULTI-TRACK™ balloon catheter, a larger rapid exchange balloon catheter or both.

Figure 1


Figure 2

Figure 3

Figure 4

Figure 5

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
## REFERENCES

1. Wilkins GT, Weyman AE, Abscal VM, Block PC, Palacios IF: Percutaneous balloon dilatation of the mitral valve: an analysis of echocardiographic variables related to outcome and the mechanism of dilatation. Br. Heart J. 60:299-308, (1988).
2. Bonhoeffer, Philipp, Piechaud JF, Sidi D, Yonga G, Jowi C, Joshi M, Mugo M, Kachaner J, Parenzan L: Mitral dilatation with the Multi-Track System: An alternative approach. Catheterization and Cardiovascular Diagnosis 36:189-193 (1995).
3. Bonhoeffer, Philipp: Percutaneous Mitral Valve Dilatation with the Multi-Track System, Catheterization and Cardiovascular Interventions, 48:178-183, (1999).
4. Tegtmeier, Charles J., M.D. & Bezirdijan Diran R., M.D. " Removing the Stuck, Ruptured Angioplasty Balloon catheter." Radiology, Vol. 139, 132-232, (April 1981).
5. Bonhoeffer, Philipp; Piechaud, Jean-Francois; The Multi-Track Angiography Catheter: A new tool for complex catheterisation in congenital heart disease, Heart, 76:173-177, (1996).

This page intentionally left blank.

This page intentionally left blank.

	Do Not Reuse Ne pas réutiliser Non riutilizzare Nicht wiederverwenden No reutilizar	Får ej återanvändas Må ikke genbruges Niet hergebruiken NÃO Reutilizar Μην επαναχρησιμοποιείτε		Catalogue Number Número de catalogue Numero di catalogo Καταλόγnummer Número de catálogo	Katalognummer Καταλόγnummer Catalogusnummer Número de Catálogo Αριθμός καταλόγου
	Batch Code Code du lot Codice del lotto Chargencode Código del lote	Sändningskod Partikode Particode Código do Lote Κωδικός παρτίδας		Use By À utiliser avant Utilizzare entro Verwendbar bis Usar antes de	Använd senast Användes inden Gebruiken voor Usar Até Ημερομηνία λήξης
	Do not resterilize Ne pas restériliser Non risterilizzare Nicht reesterilieren No reesterilizar			Får inte omsteriliseras Må ikke resteriliseres Niet hersteriliseren NÃO reesterilizar Μην επαναποστειρώνετε	
	Attention: See Instructions For Use Attention : Consulter le mode d'emploi Attenzione – Consultare le Istruzioni per l'uso Achtung: Siehe Gebrauchsanweisung Atención: consultar instrucciones de uso			OBS! Se bruksanvisningen Bemærk: Se brugsanvisningen Let op: zie gebruiksaanwijzing Atenção: Consulte as instruções de utilização Προσοχή: Δείτε τις οδηγίες χρήσης	
	Keep Away From Sunlight Garder à l'abri de la lumière du soleil Proteggere dai raggi solari Sonneneinstrahlung vermeiden Mantener alejado de la luz solar			Exponera inte för direkt solljus Holdes væk fra sollys Uit direct zonlicht houden Mantel Afstando da Luz Solar Διατηρείτε μακριά από ηλιακό φως	
	Temperatur Limitation Limite de température Limitazione di temperatura Temperaturbegrenzung Limitación de la temperatura			Temperaturbegrensning Temperaturbegrænsning Temperaturberøking Limite de Temperatura Περιορισμός θερμοκρασίας	
	Manufacturer / Date of Manufacture Fabricant / Date de fabrication Produttore / Data di fabbricazione Hersteller / Herstellungsdatum Fabricante / Fecha de fabricación			Tilværkare / Tilværkningsdato Producent / fremstillingsdato Fabrikant / Productiedatum Fabricante / Data de Fabrico Κατασκευαστής / Ημερομηνία κατασκευής	
	Sterilized Using Ethylene Oxide Stérilisé à l'oxyde d'éthylène Sterilizzato con ossido di etilene Sterilisiert mit Ethylenoxid Esterilizado con óxido de etileno			Steriliserad med etylenoxid Steriliseret med ethylenoxid Gesteriliseerd met ethyleenoxide Esterilizado Usando Óxido de Etileno Αποστειρώθηκε με χρήση οξείδιου του αιθυλενίου	
	Authorised Representative in the European Community Représentant autorisé dans la Communauté européenne Rappresentante autorizzato nella Comunità Europea Autorisierter Vertreter in der Europäischen Gemeinschaft Representante autorizado en la Comunidad Europea			Auktoriserad representant inom EU Autoriseret repræsentant i Det Europæiske Fællesskab Gevoelmsachtigd vertegenwoordiger in de Europese Gemeenschap Representante Autorizado na Comunidade Europeia Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα	
	Do not use if package is damaged Ne pas utiliser si l'emballage est endommagé Non utilizzare se l'imballaggio non è integro Nicht verwenden, wenn die Verpackung beschädigt ist No utilizar si el envase está dañado			Får inte användas om förpackningen är skadad Må ikke anvendes, hvis emballagen er beskadiget Niet gebruiken als de verpakking beschadigd is NÃO utilizar se a embalagem estiver danificada Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά	
	MR Conditional Utilisable(s) en environnement IRM sous certaines conditions A compatibilità RM condizionata Bedingt MR-sicher Producto seguro bajo ciertas condiciones de la resonancia magnética			Vilkorligt godkänd för MR Betinget MR-sikker MRI-veilig onder bepaalde voorwaarden MR Conditional (pode ser submetido a exames de RM sob determinadas condições) Ασφαλές για μαγνητική τομογραφία υπό προϋποθέσεις	
	Medical Device Dispositif Médical Dispositivo médico Medizinprodukt Producto sanitario			Medicinteknisk produkt Medicinsk udstyr Medisch hulpmiddel Dispositivo médico Ιατροτεχνολογικό προϊόν	

IFU-390CE-22

NuMED, Inc.  
2880 Main Street  
Hopkinton, NY  
USA 12965  
TEL: (315) 328-4491  
FAX: (315) 328-4941  
email: info@numedusa.com  
www.numedforchildren.com

European Representative:  
G. van Wageningen B.V.  
Hallengeweg 40, 5683 CT Best,  
The Netherlands  
TEL: +31-499-377388  
FAX: +31-499-377456  
email: info@heartmedical.nl

RM0344-181

12 July 2021