
CE
1639

PTS-X™

SIZING BALLOON CATHETER

Instructions for Use	2
Mode d'emploi.....	4
Istruzioni per l'uso	7
Gebrauchsanweisung.....	10
Instrucciones de utilización.....	13
Bruksanvisning.....	16
Brugsanvisning.....	18
Gebruiksaanwijzing	20
Instruções de utilização.....	23
Οδηγίες χρήσης.....	26

CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician.


Instructions for Use

INDICATIONS: Recommended for use in those patients with cardiovascular defects wherein accurate measurement of the defect is important to select the appropriately sized occluder device.

DESCRIPTION

The NuMED PTS-X Sizing Balloon Catheter is a coaxial design catheter with a balloon mounted on its distal tip. The lumen labeled with the balloon size is for balloon inflation while the through lumen allows the catheter to track over a guidewire. This lumen has two radiopaque platinum marker bands under the balloon shoulders and there are two additional radiopaque platinum marker bands spaced at 10 mm (as measured from leading edge to leading edge). These bands are located at the balloon center and are used as a distance reference (the 1cm balloon lengths have only two radiopaque marker bands under the balloon shoulders). Each balloon inflates to the stated diameter and length at a specific pressure. The balloon size is $\pm 10\%$ at the Rated Burst Pressure (RBP). The RBP is different for each size. Check the package label for the RBP. It is important that the balloon not be inflated beyond RBP.

**HOW SUPPLIED**

Supplied sterilized by ethylene oxide gas. Sterile and non-pyrogenic if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

WARNINGS

- CAUTION: Do not exceed the RBP. An inflation device with pressure gauge is recommended to monitor pressure. Pressure in excess of the RBP can cause balloon rupture and potential inability to withdraw the catheter through the introducer sheath.
- Use only appropriate balloon inflation medium. Do not use air or gaseous medium to inflate the balloon.
- Do not advance the guidewire, balloon dilatation catheter, or any other component if resistance is met, without first determining the cause and taking remedial action.
- This catheter is not recommended for pressure measurement or fluid injection.
- Do not remove the guidewire from the catheter at any time during the procedure.
- This device is intended for single use only. Do not resterilize and/or reuse it, as this can potentially result in compromised device performance and increased risk of cross-contamination.

PRECAUTIONS

- One should always select a PTS-X diameter larger than the unstretched defect diameter, i.e., TEE ASD size 12 mm – select 20 or 25 mm PTS-X.
- Caution should be used when inflating the balloon, over inflation can cause trauma and overstretching of the septum.
- Sizing procedures should be conducted under fluoroscopic guidance with appropriate x-ray equipment.
- Guidewires are delicate instruments. Care should be exercised while handling to help prevent the possibility of breakage.
- Careful attention must be paid to the maintenance of tight catheter connections and aspiration before proceeding to avoid air introduction into the system.
- Under no circumstances should any portion of the catheter system be advanced against resistance. The cause of the resistance should be identified with fluoroscopy and action taken to remedy the problem.
- If resistance is felt upon removal, then the balloon, guidewire, and the sheath should be removed together as a unit, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the balloon catheter and sheath as a unit and withdrawing both together, using a gentle twisting motion combined with traction.
- Before removing catheter from sheath it is very important that the balloon is completely deflated.
- Proper functioning of the catheter depends on its integrity. Care should be used when handling the catheter. Damage may result from kinking, stretching, or forceful wiping of the catheter.

POTENTIAL COMPLICATIONS/ADVERSE EFFECTS

Potential balloon separation following balloon rupture or abuse and the subsequent need to use a snare or other medical interventional techniques to retrieve the pieces.

NOTE: There have been infrequent reports of larger diameter balloons bursting circumferentially, possibly due to a combination of tight focal strictures in large vessels. In any instance of a balloon rupture while in use, it is recommended that a sheath be placed over the ruptured balloon prior to withdrawal through the entry site. This can be accomplished by cutting off the proximal end of the catheter and slipping an appropriately sized sheath over the catheter into the entry site. For specific technique, refer to: Tegtmeier, Charles J., M.D. & Bezirdjian Diran R., M.D. "Removing the Stuck, Ruptured Angioplasty Balloon Catheter." Radiology, Volume 139, 231-232, April 1981.

Potential complications & adverse effects associated with device use and indication include:

- Trauma and overstretching of the septum
- Device embolization
- Device erosion
- Access site complications

Any serious incident that has occurred in relation to the device should be reported to NuMED & the Competent Authority in the country of use.

INSPECTION AND PREPARATION

1. Insert guidewire through the distal tip until guidewire exceeds proximal port.
2. Remove balloon protector. Inspect the catheter for damage prior to insertion.
3. Perform sizing using either a 1:3 or a 1:4 solution of saline and contrast medium, respectively.
4. Attach an inflation device with pressure gauge half filled with the contrast solution to the balloon port of the catheter.
5. Purge the catheter through lumen thoroughly, observing for leaks.
6. To check inflation/deflation times, use a stopwatch. Repeat the procedure several times to verify the inflation/deflation time.
7. Point the inflation device with pressure gauge nozzle downward and aspirate until all air is removed from the balloon, and bubbles no longer appear in the contrast solution.
8. Turn the stopcock off to maintain the vacuum in the balloon.
9. Remove guidewire.

INSERTION: VASCULAR

1. Enter the vessel percutaneously using the standard Seldinger technique over a 0.035" superstiff wire positioned in the left upper pulmonary vein or left atrial appendage.
2. Position the balloon across the defect and slowly inflate with diluted contrast. Both the Transesophageal echocardiogram (TEE) and fluoroscopy monitor the inflation.
3. When the defect is completely occluded there will be no color flow shunt by TEE and a waist appears in the balloon.
4. The balloon waist size is measured by TEE or cineangiogram. This size corresponds to the stretched size of the defect.

DEFLATION AND WITHDRAWAL

1. Deflate the balloon by drawing a vacuum with an inflation device with pressure gauge. Note: The greater the vacuum applied and held during withdrawal, the lower the deflated balloon profile.
2. Gently withdraw the catheter. As the balloon exits the vessel, use a smooth, gentle, steady motion. If resistance is felt upon removal, then the balloon, guidewire, and the sheath should be removed together as a unit under fluoroscopic guidance, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the balloon catheter and sheath as a unit and withdrawing both together, using a gentle twisting motion combined with traction.
3. Apply pressure to the insertion site according to standard practice or hospital protocol for percutaneous vascular procedures.
4. Dispose of device after use according to standard hospital protocol for biohazardous devices.

WARNING: NuMED catheters are placed in the extremely hostile environment of the human body. Catheters may fail to function for a variety of causes including, but not limited to, medical complications or failure of catheters by breakage. In addition, despite the exercise of all due care in design, component selection, and testing prior to sale, catheters may be easily damaged before, during, or after insertion by improper handling or other intervening acts. Consequently, no representation or warranty is made that failure or cessation of function of catheters will not occur or that the body will not react adversely to the placement of catheters or that medical complications will not follow the use of catheters.

NuMED cannot warrant or guarantee NuMED accessories because the structure of the accessories may be damaged by improper handling before or during use. Therefore, no representations or warranties are made concerning them.

Warranty and Limitations

Catheters and accessories are sold in an 'as is' condition. The entire risk as to the quality and performance of the catheter is with the buyer. NuMED disclaims all warranties, expressed or implied, with respect to catheters and accessories, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose. NuMED shall not be liable to any person for any medical expenses or any direct or consequential damages resulting from the use of any catheter or accessory or caused by any defect, failure, or malfunction of any catheter or accessory, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind NuMED to any representation or warranty with respect to catheters and accessories.

PTS-X™ Balloon Sizing Chart

Applied Pressure	10.0 (mm)	12.0 (mm)	15.0 (mm)	20.0 (mm)	25.0 (mm)	30.0 (mm)	35.0 (mm)	40.0 (mm)
0.5 ATM	8.35	9.77	12.57	-	-	-	31.74	38.72
1.0 ATM	9.16	10.57	13.59	19.87	24.81	30.12	34.28	
1.5 ATM	9.67	11.43	14.82	20.86	26.56			
2.0 ATM	10.08	12.03	15.52					
2.5 ATM	10.37							

The figures in bold face represent the balloon diameter @ Rated Burst Pressure.

FOR ALL NuMED CATHETERS AN INFLATION DEVICE WITH PRESSURE GAUGE SHOULD BE USED.

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
REFERENCES

- Zhu W, Cao QL, Rhodes J, Hijazi ZM. Measurement of atrial septal defect size: a comparative study between three-dimensional transesophageal echocardiography and the standard balloon sizing methods. *Pediatric Cardiology* 2000 Sep-Oct;21(5):465-9.
- Acar P, Saliba Z, Bonhoeffer P, Aggoun Y, Bonnet D, Sidi D, Kachaner J. Influence of atrial septal defect anatomy in patient selection and assessment of closure with the Cardioseal device; a three-dimensional transesophageal echocardiographic reconstruction. *Eur. Heart Journal* 2000 Apr;21(7):573-81.
- GU X, Han YM, Berry J, Urness M, Amplatz K. A new technique for sizing of atrial septal defects. *Catheter Cardiovascular Intervention* 1999 jan;46(1):51-7.

This page intentionally left blank

This page intentionally left blank

	Do Not Reuse Ne pas réutiliser Non riutilizzare Nicht wiederverwenden No reutilizar	Får ej återvändas Må ikke genbruges Niet hergebruiken Não Reutilizar Μην επαναχρησιμοποιείτε		Catalogue Number Número de catalogue Numero di catalogo Katalognummer Número de catálogo	Katalognummer Catalognummer Katalogusnummer Número de Catálogo Αριθμός καταλόγου
	Batch Code Code du lot Codice del lotto Chargenode Código del lote	Sändningskod Partikode Partijcode Código do Lote Κωδικός παρτίδας		Use By À utiliser avant Utilizzare entro Verwendbar bis Usar antes de	Använd senast Används inden Gebruiken voor Usar Até Ημερομηνία λήξης
	Do not resterilize Ne pas restériliser Non resterilizzare Nicht resterilisieren No reesterilizar			Får inte omsteriliseras Må ikke resteriliseres Niet hersteriliseren Não reesterilizar Μην επανααστεριώνετε	
	Attention: See Instructions For Use Attention : Consulter le mode d'emploi Attenzione – Consultare le istruzioni per l'uso Achtung: Siehe Gebrauchsanweisung Atención: consultar instrucciones de uso			OBS! Se brugsanvisningen Bemærk: Se brugsanvisningen Let op: zie gebruiksaanwijzing Atenção: Consulte as instruções de utilização Προσοχή: Δείτε τις οδηγίες χρήσης	
	Keep Away From Sunlight Garder à l'abri de la lumière du soleil Proteggere dai raggi solari Sonneneinstrahlung vermeiden Mantener alejado de la luz solar			Exponera inte för direkt solljus Holdes væk fra sollys Uit direct zonlicht houden Manter Afastado da Luz Solar Διατηρείτε μακριά από ηλιακό φως	
	Temperature Limitation Limite de température Limitazione di temperatura Temperaturbegrenzung Limitación de la temperatura			Temperaturbegränsning Temperaturbegränsning Temperatuurbepörking Limite de Temperatura Περιορισμός θερμοκρασίας	
	Manufacturer / Date of Manufacture Fabricant / Date de fabrication Produttore / Data di fabbricazione Hersteller / Herstellungsdatum Fabricante / Fecha de fabricación			Tilverkare / Tilverkningsdatum Producent / fremstillingsdato Fabrikant / Productiedatum Fabricante / Data de Fabrico Κατασκευαστής / Ημερομηνία κατασκευής	
	Sterilized Using Ethylene Oxide Stérilisé à l'oxyde d'éthylène Sterilizzato con ossido di etilene Sterilisiert mit Ethylenoxid Esterilizado con óxido de etileno			Steriliserad med etylenoxid Steriliseret med ethylenoxid Gesteriliseerd met ethylenoxide Esterilizado Usando Óxido de Etileno Αποστειρώθηκε με χρήση οξειδίου του αιθυλενίου	
	Authorised Representative in the European Community Représentante autorisé dans la Communauté européenne Rappresentante autorizzato nella Comunità Europea Autorisierter Vertreter in der Europäischen Gemeinschaft Representante autorizado en la Comunidad Europea			Auktoriserad representant inom EU Autoriseret repræsentant i Det Europæiske Fællesskab Gevolmachtigd vertegenwoordiger in de Europese Gemeenschap Representante Autorizado na Comunidade Europeia Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα	
	Do not use if package is damaged Ne pas utiliser si l'emballage est endommagé Non utilizzare se l'imballaggio non è integro Nicht verwenden, wenn die Verpackung beschädigt ist No utilizar si el envase está dañado			Får inte användas om förpackningen är skadad Må ikke anvendes, hvis emballagen er beskadiget Niet gebruiken als de verpakking beschadigd is Não utilizar se a embalagem estiver danificada Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά	
	Medical Device Dispositif Médical Dispositivo médico Medizinprodukt Producto sanitario			Medicinteknisk produkt Medicinsk udstyr Medisch hulpmiddel Dispositivo médico Ιατροτεχνολογικό προϊόν	

IFU-360XC-13

NuMED, Inc.
2880 Main Street
Hopkinton, NY
USA 12965
TEL: (315) 328-4491
FAX: (315) 328-4941
email: info@numedusa.com
www.numedforchildren.com

European Representative:
G. van Wageningen B.V.
Hallenweg 40, 5683 CT Best,
The Netherlands
TEL: +31-499-377388
FAX: +31-499-377456
email: info@heartmedical.nl

RM0344-27L

02 August 2021