
CE  
1639

# TEARDUCT™ CATHETER

## DACRYOCYSTOPLASTY BALLOON DILATATION

Instructions for Use .....	2
Mode d'emploi.....	4
Istruzioni per l'uso .....	6
Gebrauchsanweisung.....	8
Instrucciones de utilización.....	10
Bruksanvisning.....	12
Brugsanvisning .....	14
Gebruiksaanwijzing .....	16
Instruções de utilização.....	18
Οδηγίες χρήσης .....	20

**CAUTION:** Federal (USA) Law restricts this device to sale by or on the order of a physician.

## Instructions for Use

**INDICATIONS**

Epiphora, Obstruction of the nasolacrimal duct system.

**DESCRIPTION**

The NuMED Tearduct Catheter is a coaxial design catheter with a distally mounted tipless balloon. The lumen labeled with the balloon size is for balloon inflation while the through lumen allows the catheter to track over a guidewire. A radiopaque band defines the real tip of the dilated balloon. Each balloon inflates to the stated diameter and length at a specific pressure. The balloon size is  $\pm 10\%$  at the Rated Burst Pressure (RBP). Check the package label for the RBP. It is important that the balloon not be inflated beyond the RBP.

**HOW SUPPLIED**

Supplied sterilized by ethylene oxide gas. Sterile and non-pyrogenic if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

**CONTRAINDICATIONS**

Other than standard risk associated with insertion of a catheter, there are no known contraindications for Dacryocystoplasty procedures. The patient's medical condition could affect successful use of this catheter.

**WARNINGS**

- **CAUTION:** Do not exceed the RBP. An inflation device with pressure gauge is recommended to monitor pressure. Pressure in excess of the RBP can cause balloon rupture and potential inability to withdraw the catheter through the introducer sheath.
- Catheter balloon inflation diameter must be carefully considered in selecting a particular size for any patient. The inflated balloon diameter should not be significantly greater than the tearduct. It is important to perform a dacryocystography prior to dilatation to measure the size of the tearduct.
- Use only appropriate balloon inflation medium. Do not use air or gaseous medium to inflate the balloon.
- This catheter is not recommended for pressure measurement or fluid injection.
- Do not remove the guidewire from the catheter at any time during the procedure.
- This device is intended for single use only. Do not resterilize and/or reuse it, as this can potentially result in compromised device performance and increased risk of cross-contamination.
- The catheter is intended for Dacryocystoplasty only and is not intended for angioplasty.
- **THE CATHETER IS NOT INTENDED FOR USE WITH STENTS.**
- The catheter should be used prior to the 'Use Before' date noted on the package label.

**PRECAUTIONS**

- Dilatation procedure should be conducted under fluoroscopic guidance with appropriate x-ray equipment.
- Guidewires are delicate instruments. Care should be exercised while handling to help prevent the possibility of breakage.
- Careful attention must be paid to the maintenance of tight catheter connections and aspiration before proceeding to avoid air introduction into the system.
- Under no circumstances should any portion of the catheter system be advanced against resistance. The cause of the resistance should be identified with fluoroscopy and action taken to remedy the problem.
- If resistance is felt upon removal, then the balloon and guidewire (total system) should be removed together as a unit, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the balloon catheter and guidewire as a unit and withdrawing both together, using a gentle twisting motion combined with traction.
- Before removing the catheter it is very important that the balloon is completely deflated.
- Proper functioning of the catheter depends upon its integrity. Care should be used when handling the catheter. Damage may result from kinking, stretching or forceful wiping of the catheter.

**INSTRUCTIONS FOR USE**

Prior to use, carefully examine all equipment to be used during the procedure, including the catheter, to verify proper function and that the catheter size is suitable for the specific procedure for which it is intended. Also, inflate the dilatation catheter to the appropriate RBP and deflate to verify proper function.

1. Remove balloon protector. Inspect the catheter for damage prior to insertion.
2. Check that all the connections are tight. Fill and purge the dilatation balloon. Prime and flush the distal lumen.
3. Under fluoroscopic guidance advance the guidewire antegradely into the lacrimal drainage system via the lower or upper lacrimal point. When the guidewire extends sufficiently from the nose, the catheter is passed retrogradely over the wire into the lacrimal drainage system. A radiopaque band defines the tip of the dilatation balloon.
4. The distal lumen is provided for guidewire tracking. An inflation device with pressure gauge is required to monitor inflation pressure (refer to package label for RBP)
5. Perform dilations using either a 50/50 or a 75/25 solution of saline and contrast medium, respectively. Patient monitoring is required during dilations. Balloon can be either partially or fully inflated to achieve dilatation. **DO NOT EXCEED THE RBP.**
6. Deflate the balloon by drawing a vacuum with an inflation device with pressure gauge. Note: The greater the vacuum applied and held during withdrawal, the lower the deflated balloon profile. Gently withdraw the catheter. If resistance is felt upon removal, then the balloon and the guidewire should be removed together as a unit under fluoroscopic guidance particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the balloon and guidewire as a unit and withdrawing both together, using a gentle twisting motion combined with traction.
7. Dispose of device after use according to standard hospital protocol for biohazardous devices.

**POTENTIAL COMPLICATIONS/ADVERSE EFFECTS**

Potential balloon separation following balloon rupture or abuse and the subsequent need to use a snare or other medical interventional technique to retrieve the pieces.

Potential complications & adverse effects associated with device use and indication include:

- Perforation
- Hematoma
- Damage of the nasolacrimal system
- Inflammation
- Infection
- Irritation

Any serious incident that has occurred in relation to the device should be reported to NuMED & the Competent Authority in the country of use.

**WARNING:** NuMED catheters are placed in the extremely hostile environment of the human body. Catheters may fail to function for a variety of causes including, but not limited to, medical complications or failure of catheters by breakage. In addition, despite the exercise of all due care in design, component selection, manufacture and testing prior to sale, catheters may be easily damaged before, during, or after insertion by improper handling or other intervening acts. Consequently, no representation or warranty is made that failure or cessation of function of catheters will not occur or that the body will not react adversely to the placement of catheters or that medical complications will not follow the use of catheters.

NuMED cannot warrant or guarantee NuMED accessories because the structure of the accessories may be damaged by improper handling before or during use. Therefore, no representations or warranties are made concerning them.

**Warranty and Limitations**

Catheters and accessories are sold in an 'as is' condition. The entire risk as to the quality and performance of the catheter is with the buyer. NuMED disclaims all warranties, expressed or implied, with respect to catheters and accessories, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose. NuMED shall not be liable to any person for any medical expenses or any direct or consequential damages resulting from the use of any catheter or accessory or caused by any defect, failure, or malfunction of any catheter or accessory, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind NuMED to any representation or warranty with respect to catheters and accessories.

**Tearduct™ Balloon Sizing Chart**

Applied Press.	3.0 (mm)
1.0 ATM	2.35
2.0 ATM	2.46
3.0 ATM	2.54
4.0 ATM	2.63
5.0 ATM	2.71
6.0 ATM	2.77
7.0 ATM	2.82
8.0 ATM	2.88
9.0 ATM	2.95
10.0 ATM	3.01
11.0 ATM	3.04
12.0 ATM	3.08
13.0 ATM	3.12
<b>14.0 ATM</b>	<b>3.15</b>

The figures in bold face represent the balloon diameter @ Rated Burst Pressure.

**FOR ALL NuMED CATHETERS AN INFLATION DEVICE WITH PRESSURE GAUGE SHOULD BE USED.**

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
## REFERENCES

- Janssen AG, Manson K, Krabbe GJ, Van de Veen S, Helder AK., Dacryocystoplasty: treatment of epiphora by means of balloon dilatation of the obstructed nasolacrimal duct system. *Radiology* 1994; 193: 453-456
- Emmerich KH, Busse H, Meyer-Rusenberg HW.  
Dacryocystorhinostomia externa., Technik, Indikationen, Ergebnisse. *Ophthalmologie* 1994; 91: 395-398
- Song HY, Ahn HS, Park CK, Kwon SH, Kim CS, Choi KC., Complete obstruction of the nasolacrimal system. Part I. Treatment with balloon dilatation. *Radiology* 1993; 186: 367-371
- Lee JM, Song HY, Han YM et al., Balloon dacryocystoplasty; results in the treatment of complete and partial obstructions of the nasolacrimal system. *Radiology* 1994; 192: 503-508
- Ilgit ET, Yuksel D, Unal M, Akpek S, Isik S, Hasanreisoglu B., Transluminal balloon dilatation of the lacrimal drainage system for the treatment of epiphora. *American Journal of Radiology* 1995; 165: 1517-1524
- Janssen AG, Mansour K, Bos JJ., Obstructed nasolacrimal duct system in epiphora: long-term results of dacryocystoplasty by means of balloon dilatation. *Radiology* 1997; 205: 791-796
- Berkefeld J, Krichner J, Muller HM, Fries U, Kollath J. Balloon dacryocystoplasty: indications and contraindications. *Radiology* 1997; 205: 785-790
- Janssen AG et al., Diagnoses and treatment of epiphora due to obstruction of the nasolacrimal duct system. A radiological study, Thesis February 2000

This page intentionally left blank

	Do Not Reuse Ne pas réutiliser Non riutilizzare Nicht wiederverwenden No reutilizar Får ej återvändas	Μά ikke genbruges Niet hergebruken Não Reutilizar Tekrar Kullanmayın Μην επαναχρησιμοποιείτε		Catalogue Number Número de catalogue Número di catalogo Katalognummer Número de catálogo Katalognummer	Katalognummer Catalogusnummer Número de Catálogo Katalog Numarası Αριθμός καταλόγου
	Batch Code Code du lot Codice del lotto Chargecode Código del lote Sändningskod	Partikode Partijcode Código do Lote Parti Kodu Κωδικός παρτίδας		Use By À utiliser avant Utilizzare entro Verwendbar bis Usar antes de Använd senast	Avendes inden Gebruiken voor Usar Até Son Kullanma Tarihi Ημερομηνία λήξης
	Do not sterilize Ne pas restériliser Non sterilizzare Nicht sterilisieren No reesterilizar Får inte omsteriliseras	Μά ikke reesteriliseres Niet hersteriliseren NÃO reesterilizar Yeniden sterilize etmeyin Μην επαναστεριώνετε		Temperature Limitation Limite de température Limitazione di temperatura Temperaturbegrenzung Limitación de la temperatura Temperaturbegrensning	Temperaturbegrensning Temperaturbeperking Limite de Temperatura Isi Sınırlaması Περιορισμός θερμοκρασίας
	Consult Instructions for Use Consulter le mode d'emploi Consultare le istruzioni per l'uso Gebrauchsanweisung beachten Consultar las instrucciones de uso Läs bruksanvisningen			Konsultér brugsanvisningen Gebrauchsanweisung raadplegen Consultar as instruções de utilização Kullanım talimatına bakın Συμβουλευτείτε τις οδηγίες χρήσης	
	Keep Away From Sunlight Garder à l'abri de la lumière du soleil Proteggere dai raggi solari Sonneneinstrahlung vermeiden Mantener alejado de la luz solar Exponera inte för direkt solljus			Holdes væk fra sollys Uit direct zonlicht houden Mantar Alastado da Luz Solar Güneş Işığından Uzak Tutun Διατηρείτε μακριά από ηλιακό φως	
	Manufacturer / Date of Manufacture Fabricant / Date de fabrication Produttore / Data di fabbricazione Hersteller / Herstellungsdatum Fabricante / Fecha de fabricación Tilvarekare / Tilverkningsdatum			Product / fremstillingsdato Fabrikant / Productiedatum Fabricante / Data de Fabrico Üretici / Üretim Tarihi Κατασκευαστής / Ημερομηνία κατασκευής	
	Double Sterile Barrier System Système de double barrière stérile Sistema di barriera sterile doppia Doppeltes Sterilbarriersystem Sistema de doble barrera estéril Dubbelt sterilt barrieresystem			Dobbelt sterilt barrieresystem System met dubbele steriele barrière Sistema de barreira dupla estéril Çift steril bariyer sistemi Διπλά αποστειρωμένο σύστημα φραγμού	
	Sterilized Using Ethylene Oxide Stérilisé à l'oxyde d'éthylène Sterilizzato con ossido di etilene Sterilisiert mit Ethylenoxid Esterilizado con óxido de etilén Steriliserad med etylenoxid			Steriliseret med ethylenoxid Gesteriliseerd met ethylenoxide Esterilizado Usando Óxido de Etileno Etilen Oksit Kullanılarak Sterilize Edilmiştir Αποστειρώθηκε με χρήση οξείδιου του αιθυλενίου	
	Authorised Representative in the European Community Représentante autorisé dans la Communauté européenne Rappresentante autorizzato nella Comunità Europea Autorisierter Vertreter in der Europäischen Gemeinschaft Representante autorizado en la Comunidad Europea Auktoriserad representant inom EU			Autoriseret representant i Det Europæiske Fællesskab Gevnoamtlichgd vertegenwoordiger in de Europese Gemeenschap Representante Autorizado na Comunidade Europeia Avruppa Topluluğunda Yetkili Temsilci Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα	
	Do not use if package is damaged and consult instructions for use Ne pas utiliser si l'emballage est endommagé et consulter le mode d'emploi Non usare se la confezione è danneggiata e consultare le istruzioni per l'uso Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten No utilizar si el envase está dañado y consultar las instrucciones de uso Får ej användas om förpackningen är skadad och se bruksanvisningen			Μά ikke anvendes, hvis emballagen er beskadiget. Konsultér brugsanvisningen Niet gebruiken als de verpakking is beschadigd en gebruiksaanwijzing raadplegen Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização Paket hasarlıysa kullanmayın ve kullanım talimatına bakın Μην χρησιμοποιείτε, ένα ή συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	
	Medical Device Dispositif Médical Dispositivo medico Medizinprodukt Producto sanitario Medicinteknisk produkt			Medicinsk udstyr Medisch hulpmiddel Dispositivo médico Tibbi Cihaz Ιατροτεχνολογικό προϊόν	

IFU-352CE-11

NuMED, Inc.  
2880 Main Street  
Hopkinton, NY  
USA 12965  
TEL: (315) 328-4491  
FAX: (315) 328-4941  
email: [info@numedusa.com](mailto:info@numedusa.com)  
[www.numedforchildren.com](http://www.numedforchildren.com)

European Representative:  
G. van Wageningen B.V.  
Hallenweg 40, 5683 CT Best,  
The Netherlands  
TEL: +31-499-377388  
FAX: +31-499-377456  
email: [info@heartmedical.nl](mailto:info@heartmedical.nl)

RM0344-161

03 February 2022