
# COVERED MOUNTED CP STENT™

## COARCTATION OF THE AORTA AND RIGHT VENTRICULAR OUTFLOW TRACT

|                                    |    |
|------------------------------------|----|
| Instructions for use .....         | 2  |
| Mode d'emploi .....                | 5  |
| Istruzioni per l'uso .....         | 9  |
| Gebrauchsanweisung .....           | 13 |
| Instrucciones de utilización ..... | 17 |
| Bruksanvisning .....               | 21 |
| Brugsanvisning .....               | 25 |
| Gebruiksaanwijzing .....           | 29 |
| Instruções de utilização .....     | 33 |
| Kullanım Talimatları .....         | 37 |
| Οδηγίες χρήσης .....               | 41 |

**CAUTION:** Federal (USA) Law restricts this device to sale by or on the order of a physician.

The ACC/AHA 2008 Guidelines for the Management of Adults With Congenital Heart disease recommends a yearly follow-up and additional imaging of the Coarctation site by CT or MRI every 5 years or less.

## Instructions for Use

**INTENDED PURPOSE**

The stent is intended to dilate aortic coarctations, and/or contain conduit disruptions in the right ventricular outflow tract (RVOT), using percutaneous implantation techniques.

**INDICATIONS**Coarctation of the Aorta (CoA)

Indicated for treatment of native and/or recurrent coarctation of the aorta in patients with the following clinical conditions:

- Stenosis of the aorta resulting in significant anatomic narrowing as determined by angiography or non-invasive imaging, i.e. echocardiography, magnetic resonance imaging (MRI), CT scan;
- Stenosis of the aorta resulting in hemodynamic alterations, resulting in systolic pressure gradient, systemic hypertension or altered left ventricular function;
- Stenosis of the aorta where balloon angioplasty is ineffective or contraindicated;
- Stenosis diameter <20% of the adjacent vessel diameter. Stenosis that would present increased risk of vascular damage or disruption; or aneurysm associated with coarctation of the aorta.

Right Ventricular Outflow Tract (RVOT)

Indicated for treatment of right ventricle (right ventricular outflow tract) conduit disruptions that are identified during conduit pre-dilatation procedures performed in preparation for transcatheter pulmonary valve replacement.

**DESCRIPTION**

The covered stent is balloon expandable and intended for permanent implant, to increase the minimal luminal / coarctation diameter, reduce systolic blood pressure gradient, decrease hypertension, decrease antihypertensive medication use, and to reduce / prevent aortic wall injury in CoA applications, and to increase the rate of procedural success of TPVR in RVOT applications.

The covered stent is composed of heat treated 90% platinum / 10% iridium wire that is arranged in a "zig" pattern, laser welded at each joint, and over brazed with 24K gold. The number of zigs in a row can be varied and will impact the strength of the stent as well as the eventual expanded diameter and percent stent shortening, while the number of rows will determine the unexpanded length of the stent. The covered stent has an ePTFE covering attached to the stent framework with a cyanoacrylate adhesive. This covering acts as a fluid barrier creating a fluid tight conduit through the stent length.

The NuMED BIB® (Balloon in Balloon) Catheter is a triaxial design catheter. Two lumens are used to inflate the balloons while one lumen is for tracking over a guidewire. The radiopaque platinum marker(s) are placed beneath the 'working area' of the balloon. The inner balloon is ½ of the outer balloon diameter and 1 cm shorter. Each balloon inflates to the stated diameter and length at a specific pressure. The balloon size is ±10% at rated burst pressure (RBP). The RBP is different for each size. Check the package label for the RBP. It is important that the balloon not be inflated beyond the RBP.

The purpose of the double balloon catheter is to apply an incremental inflation for the purpose of opening a vascular channel using a balloon expandable intravascular stent. The inner balloon provides initial expansion of the stent and also acts as a tool to hold the stent in place while the outer balloon is inflated. The outer balloon is then inflated securing the stent against the vessel wall.

**IMPLANT MATERIALS**

| Component | Material               | Mass <sup>1</sup> |
|-----------|------------------------|-------------------|
| Frame     | Platinum Iridium 90/10 | 2.260g            |
|           | 24K Gold (brazing)     | 0.220g            |
| Covering  | ePTFE                  | 0.125g            |
| Adhesive  | E-5 Cyanoacrylate      | 0.005g            |

<sup>1</sup> The calculations of mass are based on the materials in the largest stent size (CVRDCP10Z60).

**HOW SUPPLIED**

Supplied sterilized by ethylene oxide gas. Sterile and non-pyrogenic if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

**CONTRAINDICATIONS – CoA and RVOT**

- Patients too small to allow safe delivery of the stent without compromise to the systemic artery used for delivery;
- Unfavorable aortic anatomy that does not dilate with high pressure balloon angioplasty (CoA only);
- Occlusion or obstruction of systemic artery precluding delivery of the stent (CoA only);
- Clinical or biological signs of infection;
- Active endocarditis;
- Known allergy to aspirin, other antiplatelet agents, or heparin (CoA only);
- Pregnancy.

**CP STENT WARNINGS**

- Radiofrequency heating during MRI scans on overlapped, 10 zig CP Stents has not been evaluated, and is not recommended.
- During the RVOT clinical study, the Medtronic Melody Valve was used for valve restoration.
- As with any type of implant, infection secondary to contamination of the stent may lead to aortitis, or abscess.
- The platinum/iridium stent may migrate from the site of implant.
- Over-stretching of the artery may result in rupture or aneurysm formation.
- The inflated diameter of the stent should at least equal the diameter of the intended implant site.
- Excessive force while crimping may weaken welds of the stent.
- Crimping the 8 zig stent on a balloon catheter smaller than 12mm, and the 10 zig on a balloon catheter smaller than 26mm, may cause damage to the stent.
- Excessive handling and manipulation of the covering while crimping the stent may cause the covering to tear off of the stent.
- Crimping the device in the opposite direction of the folds in the covering may cause the covering to catch while inserting into the hemostasis valve tool and introducer. This could cause the covering to tear off of the stent.
- Pulling the Covered stent back through the introducer and/or hemostasis valve may cause the covering to catch and tear off of the stent.

**BIB STENT PLACEMENT WARNINGS**

- Do not exceed the RBP. An inflation device with pressure gauge is recommended to monitor pressure. Pressure in excess of the

RBP can cause balloon rupture and potential inability to withdraw the catheter through the introducer sheath.

- Confirm that the distal end of the introducing sheath is at least 2.5cm back from the most proximal image markers before inflating the outer balloon. Failure to do so may stretch the outer tubing and severely hinder balloon deflation.
- Use two appropriate size inflation devices with pressure gauges for inflation.
- Do not advance the guidewire, balloon dilatation catheter, or any other component if resistance is met, without first determining the cause and taking remedial action.
- This catheter is not recommended for pressure measurement or fluid injection.
- Do not remove the guidewire from the catheter at any time during the procedure except when the procedure has been completed.
- This device is intended for single use only. Do not resterilize and/or reuse it, as this can potentially result in compromised device performance and increased risk of cross contamination.

#### PRECAUTIONS

- Use of an inflation device with pressure gauge is highly recommended during this procedure.
- Stents are delicate devices. Care should be exercised while handling to help prevent the possibility of breakage.
- The stent is rigid and may make negotiation through vessels difficult.
- Dilatation procedures should be conducted under fluoroscopic guidance with appropriate x-ray equipment.
- Guidewires are delicate instruments. Care should be exercised while handling to help prevent the possibility of breakage.
- Careful attention must be paid to the maintenance of tight catheter connections and by aspiration before proceeding to avoid air introduction into the system.
- The inflation diameter of the balloon used during stent delivery should approximate the diameter of the obstructive vessel and the intended implant site.
- Under no circumstances should any portion of the catheter system be advanced against resistance. The cause of the resistance should be identified with fluoroscopy and action taken to remedy the problem.
- If resistance is felt upon removal, then the balloon, guidewire and the sheath should be removed together as a unit, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the balloon catheter and sheath as a unit and withdrawing both together, using a gentle twisting motion combined with traction.
- The balloons must be completely deflated before retracting into the sheath.
- Proper functioning of the catheter depends on its integrity. Care should be used when handling the catheter. Damage may result from kinking, stretching, or forceful wiping of the catheter.

#### POTENTIAL COMPLICATIONS/ADVERSE EFFECTS

**NOTE:** Circumferential tear of the delivery balloon catheter prior to complete expansion of the stent may cause the balloon to become tethered to the stent, requiring surgical removal. In case of rupture of an adequately sized balloon after stent expansion, it can be withdrawn and a new balloon catheter exchanged over a guidewire to complete expansion of the stent.

Cardiac catheterization carries certain risks. Potential complications & adverse effects associated with device use and indication include:

- |                         |                |                             |                            |
|-------------------------|----------------|-----------------------------|----------------------------|
| • Femoral Artery Injury | • Hematoma     | • Endocarditis              | • Sepsis / Infection       |
| • Stent Migration       | • Thrombosis   | • Stent Stenosis            | • Transitory Arrhythmia    |
| • Stent Fracture        | • Embolization | • Aneurysm / Pseudoaneurysm | • Bleeding                 |
| • Aortic Rupture / Tear | • Death        | • Stent Malposition         | • Cerebrovascular Incident |

Any serious incident that has occurred in relation to the device should be reported to NuMED & the Competent Authority in the country of use.

#### MRI SAFETY INFORMATION

Nonclinical testing and modeling has demonstrated that the CP Stent is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3 T
- Maximum spatial gradient magnetic field of 2500 gauss/cm (25 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning (Normal Operating Mode)

Based on nonclinical testing and modeling, under the scan conditions defined above, the CP Stent is expected to produce a maximum in vivo temperature rise of less than 2°C after 15 minutes of continuous scanning.

MR image quality may be compromised if the area of interest is in the same area, or relatively close to the position of the device. In nonclinical testing, the image artifact caused by the device extends approximately 3 mm from the CP Stent when imaged with a spin echo pulse sequence and 6 mm when imaged with a gradient echo pulse sequence and a 3 T MRI System. The lumen of the device was obscured.

The presence of other implants or medical circumstances of the patient may require lower limits on some or all of the above parameters.

- (1) 10 zig, 6.0 cm stent was evaluated for Magnetic Field Interaction and Image Artifact.
- (2) 8 zig, 4.5 cm stents with 10% overlap, expanded to 12 mm, were evaluated for radiofrequency heating during MRI scans. Radiofrequency heating during MRI scans on overlapped, 10 zig CP Stents has not been evaluated, and is not recommended.

#### INSTRUCTIONS FOR USE

##### Select Stent Size

1. Measure the length of the target stricture to determine the length of stent required. Size the stent length to extend slightly proximal and distal to the stricture.
2. The appropriate stent length should be selected based on covering the entire obstructed segment with a single stent.  
Note: Should more than one stent be required, place the stent most distal from the puncture site first, followed by placement of the proximal stent in tandem.
3. Measure the diameter of the reference stricture and vessel proximal and distal to the target lesion to determine the appropriate size stent and delivery system.

**NOTE:** Average stent recoil ranges from 2.1% – 4.9%.

##### Preparation of Stent Delivery System

- Visually inspect the balloon/stent assembly to assure proper placement of the stent and to assure that the stent is evenly crimped. Uneven crimping could cause the stent to deploy in a non-symmetrical manner (MAY BE CONFIRMED BY FLUOROSCOPY).

##### Stent Deployment

1. Sterile preparation and draping of the access site should be performed according to standard practice or hospital protocol.

- Use of the tools supplied with the stent is necessary to defeat the hemostasis valve without damaging the stent or covering. Once the stent is past the hemostasis valve, the tool must be pulled out of the valve.**
- The assembly is advanced through the long delivery sheath and over the stiff guidewire into the desired location for implant.
- After correct positioning of the stent, pull back on the sheath to expose the stent. Confirm proper stent position by a small injection of contrast through the sidearm of sheath or through a second catheter.
- Expand the stent initially by inflating the inner balloon until it is fully expanded. One may "reposition" the stent at this point by moving the BIB® catheter. The unexpanded outer balloon and expanded inner balloon hold the stent tightly against the BIB® catheter. DO NOT deflate the inner balloon before expansion of the outer balloon. This could cause the stent to slip off the balloon catheter.
- Confirm positioning and inflate the outer balloon to rated diameter. Do not exceed the manufacturer's balloon rated burst pressure.

#### **Delivery System Withdrawal**

- Once the stent is expanded, deflate both balloons completely and rotate to ensure the stent is free and properly deployed. If there is a residual waist in the stent, expand only the outer balloon again, making sure not to exceed the rated burst pressure.
- Remove the balloon catheter and confirm the result with angiography.
- Apply pressure to the insertion site according to standard practice or hospital protocol.
- Dispose of device after use according to standard hospital protocol for biohazardous devices.

**NOTE:** Diameter of the stent may be increased after placement by expanding with a larger diameter balloon. Do not exceed the maximum recommended expanded stent diameter of 24mm for the 8 zig stents, and 30mm for the 10 zig stents.

#### **RETURN OF EXPLANTED DEVICE**

NuMED, Inc. is interested in obtaining recovered CP Stents. Place the explanted device in a container or vial immediately after excision. For further instructions on the return of an explanted device, contact RA Manager, NuMED, Inc. 2880 Main Street, Hopkinton, New York, 12965. Phone number: 315-328-4491.

**WARNING:** NuMED stents are placed in the extremely hostile environment of the human body. Stents may fail to function for a variety of causes including, but not limited to, medical complications or failure of stent by fracture and embolization. In addition, despite the exercise of all due care in design, component selection, manufacture, and testing prior to sale, stents may be easily damaged before, during, or after insertion by improper handling, crimping or other intervening acts.

#### **Warranty and Limitations**

Stents and accessories are sold in an 'as is' condition. The entire risk as to the quality and performance of the stent is with the buyer. NuMED disclaims all warranties, expressed or implied, with respect to catheters and accessories, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose. NuMED shall not be liable to any person for any medical expenses or any direct or consequential damages resulting from the use of any catheter or accessory or caused by any defect, failure, or malfunction of any catheter or accessory, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind NuMED to any representation or warranty with respect to catheters and accessories.

**CP Stent™ Foreshortening Chart**

| Inflated Balloon Diameter | CP8Z16<br>(Stent length after expansion)<br>% Shortening | CP8Z22<br>(Stent length after expansion)<br>% Shortening | CP8Z28<br>(Stent length after expansion)<br>% Shortening | CP8Z34<br>(Stent length after expansion)<br>% Shortening | CP8Z39<br>(Stent length after expansion)<br>% Shortening |
|---------------------------|----------------------------------------------------------|----------------------------------------------------------|----------------------------------------------------------|----------------------------------------------------------|----------------------------------------------------------|
| 12mm                      | (1.61) cm<br>2.8%                                        | (2.18) cm<br>0.8%                                        | (2.62) cm<br>4.4%                                        | (3.23) cm<br>3.1%                                        | (3.72) cm<br>1.9%                                        |
| 14mm                      | (1.54) cm<br>6.5%                                        | (2.08) cm<br>5.4%                                        | (2.56) cm<br>6.8%                                        | (3.15) cm<br>5.4%                                        | (3.66) cm<br>3.6%                                        |
| 15mm                      | (1.51) cm<br>8.5%                                        | (2.02) cm<br>7.9%                                        | (2.51) cm<br>8.6%                                        | (3.10) cm<br>7.0%                                        | (3.54) cm<br>6.6%                                        |
| 16mm                      | (1.48) cm<br>10.6%                                       | (1.98) cm<br>10.1%                                       | (2.45) cm<br>10.7%                                       | (3.00) cm<br>9.8%                                        | (3.48) cm<br>8.2%                                        |
| 18mm                      | (1.43) cm<br>13.7%                                       | (1.89) cm<br>14.0%                                       | (2.38) cm<br>13.3%                                       | (2.88) cm<br>13.5%                                       | (3.20) cm<br>15.6%                                       |
| 20mm                      | (1.32) cm<br>20.0%                                       | (1.80) cm<br>17.9%                                       | (2.30) cm<br>16.3%                                       | (2.63) cm<br>20.9%                                       | (2.96) cm<br>21.9%                                       |
| 22mm                      | (1.23) cm<br>25.4%                                       | (1.67) cm<br>23.9%                                       | (2.09) cm<br>24.0%                                       | (2.46) cm<br>26.0%                                       | (2.85) cm<br>25.0%                                       |
| 24mm                      | (1.05) cm<br>36.4%                                       | (1.46) cm<br>33.8%                                       | (1.91) cm<br>30.3%                                       | (2.07) cm<br>37.9%                                       | (2.27) cm<br>40.1%                                       |

| Inflated Balloon Diameter | CP8Z45<br>(Stent length after expansion)<br>% Shortening | CP8Z50<br>(Stent length after expansion)<br>% Shortening | CP8Z55<br>(Stent length after expansion)<br>% Shortening | CP8Z60<br>(Stent length after expansion)<br>% Shortening |
|---------------------------|----------------------------------------------------------|----------------------------------------------------------|----------------------------------------------------------|----------------------------------------------------------|
| 12mm                      | (4.17) cm<br>3.8%                                        | (4.71) cm<br>6.2%                                        | (5.25) cm<br>5.0%                                        | (5.84) cm<br>4.5%                                        |
| 14mm                      | (3.97) cm<br>8.4%                                        | (4.58) cm<br>8.7%                                        | (5.11) cm<br>7.6%                                        | (5.67) cm<br>7.3%                                        |
| 15mm                      | (3.94) cm<br>9.2%                                        | (4.50) cm<br>10.3%                                       | (4.98) cm<br>10.0%                                       | (5.55) cm<br>9.2%                                        |
| 16mm                      | (3.84) cm<br>11.4%                                       | (4.42) cm<br>11.9%                                       | (4.91) cm<br>11.2%                                       | (5.43) cm<br>11.2%                                       |
| 18mm                      | (3.71) cm<br>14.5%                                       | (4.21) cm<br>16.1%                                       | (4.70) cm<br>15.1%                                       | (5.20) cm<br>14.9%                                       |
| 20mm                      | (3.27) cm<br>24.7%                                       | (3.96) cm<br>21.0%                                       | (4.43) cm<br>20.0%                                       | (4.92) cm<br>19.5%                                       |
| 22mm                      | (3.15) cm<br>27.3%                                       | (3.71) cm<br>26.0%                                       | (4.09) cm<br>26.1%                                       | (4.55) cm<br>25.5%                                       |
| 24mm                      | (2.83) cm<br>34.9%                                       | (3.33) cm<br>33.5%                                       | (3.72) cm<br>32.8%                                       | (4.14) cm<br>32.3%                                       |

| Inflated Balloon Diameter | CP10Z39<br>(Stent length after expansion)<br>% Shortening | CP10Z45<br>(Stent length after expansion)<br>% Shortening | CP10Z50<br>(Stent length after expansion)<br>% Shortening | CP10Z55<br>(Stent length after expansion)<br>% Shortening | CP10Z60<br>(Stent length after expansion)<br>% Shortening |
|---------------------------|-----------------------------------------------------------|-----------------------------------------------------------|-----------------------------------------------------------|-----------------------------------------------------------|-----------------------------------------------------------|
| 26mm                      | (3.17) cm<br><b>18.33%</b>                                | (3.44) cm<br><b>22.09%</b>                                | (4.10) cm<br><b>17.34%</b>                                | (4.24) cm<br><b>23.32%</b>                                | (4.85) cm<br><b>20.20%</b>                                |
| 28mm                      | (2.96) cm<br><b>23.68%</b>                                | (3.24) cm<br><b>26.75%</b>                                | (3.71) cm<br><b>25.11%</b>                                | (4.00) cm<br><b>27.58%</b>                                | (4.39) cm<br><b>27.87%</b>                                |
| 30mm                      | (2.58) cm<br><b>33.45%</b>                                | (3.09) cm<br><b>30.16%</b>                                | (3.26) cm<br><b>34.34%</b>                                | (3.64) cm<br><b>34.17%</b>                                | (4.11) cm<br><b>32.55%</b>                                |

CP Stent™ 8 Zig Balloon Sizing Chart

| Inner Balloon Pressure (atm) | Stent ID (mm)           |                         |                         |                         |                         |                         |                         |                         |
|------------------------------|-------------------------|-------------------------|-------------------------|-------------------------|-------------------------|-------------------------|-------------------------|-------------------------|
|                              | 12mm Diameter RBP = 7.0 | 14mm Diameter RBP = 6.0 | 15mm Diameter RBP = 5.0 | 16mm Diameter RBP = 5.0 | 18mm Diameter RBP = 4.0 | 20mm Diameter RBP = 4.0 | 22mm Diameter RBP = 3.0 | 24mm Diameter RBP = 3.0 |
| 1.0                          | 2.75                    | 3.22                    | 3.49                    | 3.75                    | 3.94                    | 4.02                    | 4.20                    | 4.28                    |
| 2.0                          | 2.85                    | 3.32                    | 3.59                    | 3.85                    | 4.36                    | 4.13                    | 4.33                    | 4.50                    |
| 3.0                          | 5.85                    | 6.91                    | 6.89                    | 7.79                    | 8.54                    | 9.20                    | 10.16                   | 10.57                   |
| 4.0                          | 6.12                    | 7.00                    | 7.02                    | 7.95                    | 8.71                    | 9.63                    | 10.40                   | 11.08                   |
| 4.5                          |                         |                         |                         |                         |                         |                         | <b>10.84</b>            | <b>11.94</b>            |
| 5.0                          | <b>6.20</b>             | <b>7.08</b>             | <b>7.10</b>             | <b>8.04</b>             | <b>8.91</b>             | <b>10.00</b>            |                         |                         |
| Outer Balloon Pressure (atm) |                         |                         |                         |                         |                         |                         |                         |                         |
| 1.0                          | 10.73                   | 13.08                   | 13.45                   | 14.87                   | 16.85                   | 17.91                   | 20.52                   | 22.79                   |
| 2.0                          | 10.86                   | 13.27                   | 14.16                   | 15.10                   | 17.06                   | 18.38                   | 21.46                   | 23.95                   |
| 3.0                          | 11.15                   | 13.50                   | 14.55                   | 15.68                   | 17.64                   | 19.42                   | <b>21.98</b>            | <b>24.68</b>            |
| 4.0                          | 11.33                   | 13.68                   | 14.88                   | 15.93                   | <b>18.06</b>            | <b>20.07</b>            |                         |                         |
| 5.0                          | 11.62                   | 13.87                   | <b>15.06</b>            | <b>16.19</b>            |                         |                         |                         |                         |
| 6.0                          | 11.80                   | <b>13.98</b>            |                         |                         |                         |                         |                         |                         |
| 7.0                          | <b>12.04</b>            |                         |                         |                         |                         |                         |                         |                         |

CP Stent™ 10 Zig Balloon Sizing Chart

| Inner Balloon Pressure (atm) | Stent ID (mm)           |                         |                         |
|------------------------------|-------------------------|-------------------------|-------------------------|
|                              | 26mm Diameter RBP = 3.0 | 28mm Diameter RBP = 2.0 | 30mm Diameter RBP = 2.0 |
| 1.0                          | 10.25                   | 10.94                   | 11.96                   |
| 2.0                          | 10.77                   | 11.39                   | 12.42                   |
| 3.0                          | 11.27                   | 11.87                   | 12.89                   |
| 4.0                          | <b>12.05</b>            | <b>12.97</b>            | <b>13.81</b>            |
| Outer Balloon Pressure (atm) |                         |                         |                         |
| 0.5                          |                         | 22.85                   | 24.84                   |
| 1.0                          | 21.62                   | 23.87                   | 25.80                   |
| 1.5                          |                         | 24.87                   | 26.81                   |
| 2.0                          | 23.34                   | <b>27.44</b>            | <b>29.94</b>            |
| 3.0                          | <b>25.44</b>            |                         |                         |

CP Stent™

| BIB DELIVERY CATHETER BALLOON DIAMETER AND INTRODUCER SIZE | REQUIRED INTRODUCER WITH COVERED CP STENT |
|------------------------------------------------------------|-------------------------------------------|
| 12MM – 14MM (8F)                                           | 12F                                       |
| 15MM – 16MM (9F)                                           | 12F                                       |
| 18MM – 20MM (10F)                                          | 14F                                       |
| 22MM – 24MM (11F)                                          | 14F                                       |
| 26MM (16F)                                                 | 16F                                       |
| 28MM – 30MM (16F)                                          | 18F                                       |

\*This data is based on testing performed using the NuMED BIB® Stent Placement Catheter.  
The figures in bold face represent the stent ID @ Rated Burst Pressure.

FOR ALL NuMED CATHETERS AN INFLATION DEVICE WITH PRESSURE GAUGE SHOULD BE USED.

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